




BD[®] OMICS-One XT WTA Assay

Automated on the Hamilton[™] Microlab[™] NGS STAR Assay Ready Workstation



Contents

-  Overview
-  Workflow and automated process
-  System description
-  Chemistry
-  Application performance
-  BOM and pricing
-  Support

Overview

Overview

Single-cell analysis

Obtaining high-quality libraries from single cells to evaluate transcriptomic heterogeneity at the single-cell level is complicated for many reasons, including:



Manual library prep is time consuming

Hours spent on library prep take away from data analysis and interpretation



Manual prep is error prone

Variability and mistakes can occur at any step, impacting data quality



Scaling is difficult


Manual methods struggle to keep pace with increasing sequencing demands

The automated BD[®] OMICS-One XT WTA Assay minimizes hands-on time and technical variability caused by human error, increases throughput to keep pace with sequencing demands and ensures consistent and quality libraries using robust reagents and pre-optimized methods.

Workflow and automated process

Workflow


BD Rhapsody™ Single-Cell Analysis System workflow




Workflow

WTA library prep workflow on the Hamilton™ Microlab™ NGS STAR Assay Ready Workstation (ARW)

Random Priming and Extension (RPE) and RPE PCR



WTA Index PCR



 Safe Stopping Point

 Automated Step

 Manual Step


System description

Features

WTA library prep from single cells on the Hamilton™ Microlab™ NGS STAR Assay Ready Workstation (ARW)



Biologically validated methods and kits



Prepare up to 48 sequencing ready libraries in a single run




Minimize human error and variability to ensure consistent, quality libraries



Pre-optimized, standardized methods with safe stopping points




Easily integrates with existing NGS STAR ARWs



Deck layout

WTA library prep from single cells on the Hamilton™ Microlab™ NGS STAR Assay Ready Workstation (ARW)



- 1 96-Well On-Deck Thermal Cycler*
- 2 Stacked Plate Lids
- 3 Stacked Plates
- 4 50-µL and 1,000-µL Conductive Tips
- 5 Tip Adapter*
- 6 300-µL Conductive Tips
- 7 60-mL Reagent Troughs
- 8 Microtubes
- 9 Plate Positions
- 10 Plate Magnet
- 11 Heater Shakers
- 12 Cold Plate (Air Cooled)
- 13 Cold Plate (Air Cooled)
- 14 Heater Shaker
- 15 Cold Plate (Air Cooled)
- 16 CO-RE® Paddles
- 17 Channel Waste
- 18 Autoload


*On-Deck Thermal Cycling and 96 Multi-Probe Head are optional

Timing and consumables requirements

WTA library prep from single cells on the Hamilton™ Microlab™ NGS STAR ARW

This solution can process 24 libraries from single cell cDNA in under 2 days

Labware	BD® OMICS-One XT WTA Assay Kit		
	8 samples/run	24 samples/run	48 samples/run
60-mL reagent reservoir	6	6	6
50-µL filtered tips	120	360	720
300-µL filtered tips	224	672	1344
1,000-µL filtered tips	192	528	1032
96-well PCR plate	4	4	4
96-well midi plate	5	5	5
2-mL screw cap tubes	Up to 23	Up to 23	Up to 23
Comfort lid	2	2	2
*Estimated hands on time	30 min	30 min	30 min
Estimated instrument processing time	7 h 17 min	9 h 17 min	14 h 50 min
Estimated total workflow time	7 h 47 min	9 h 47 min	15 h 20 min



Note: Max throughput is 48 samples

*Hands on time includes loading the instrument with consumables and reagents, not manual reagent prep time

Chemistry

Chemistry

BD® OMICS-One XT WTA Assay

Streamline your library prep workflow with the BD® OMICS-One XT WTA Assay and make automation set up easier.

- Optimized and consolidated reagents for single-cell library preparation
- Tailored for compatibility with liquid handling automation systems
- Modular workflow with safe stopping points
- Compatible with Illumina™ and Element™ sequencing platforms



Chemistry

BD® OMICS-One XT WTA Assay Kits

BD® OMICS-One XT WTA Products	Cat. No.	List Price (USD)
BD® OMICS-One XT WTA Amplification Kit (48 Rxns)	667246	\$11,502
BD® OMICS-One XT Dual Index Kit A (48 Rxns)	571973	\$450
BD® OMICS-One XT Magnetic Baseplate	667393	\$1,400

Companion Products	Cat. No.	List Price (USD)
BD Rhapsody™ HT Xpress Package	666625	\$15,000
BD Rhapsody™ Scanner	633701	\$54,450
BD Rhapsody™ 8-Lane Cartridge	666262	\$2,000
BD Rhapsody™ Enhanced Cartridge Reagent Kit V2	664887	\$3,210
BD Rhapsody™ Enhanced Cartridge Reagent Kit V3	667052	\$3,210
BD Rhapsody™ cDNA Kit	633773	\$1,550



Application performance

WTA NGS data from single cells in PBMC samples

Library DNA yield and QC metrics





Figure legend: Three runs of 18 libraries (six each of 6K subsampled beads, 6K full sampled beads and 20K full sampled beads from the same pool of cellular transcripts captured on the BD Rhapsody™ Enhanced Cell Capture Beads) were prepared on the Hamilton™ Microlab™ NGS STAR ARW along with five manual bench controls of each condition. Shown are sequenced libraries of 15 libraries total for 6K sub and 6K full each and nine libraries total for 20K full sampled beads.

A) Final quality control of indexed libraries show automated libraries was comparable to manual bench controls for Qubit concentration (ng/μL). **B)** Final quality control by TapeStation analysis shows similar average bp size distribution between 200–1,000bp.

WTA NGS data from single cells in PBMC samples


Final library QC TapeStation traces

6K sub




4 [File 1] LP081 Run2 Index Deck - 2024-06-11 - 11.53.54.HSD5000
A4 — A4_6ksub
E4 — F4_6ksub
G4 — A5_6ksub
D5 — G5_6k sub
E5 — A6_6k sub
F5 — C6_6k sub
4 [File 2] LP081 Run2 Index Bench Ctrls - 2024-06-10 - 16.19.51.HSD5000
B1 — 1_Ctrl_6ksub
C1 — 2_Ctrl_6ksub

6K full



4 [File 1] LP081 Run2 Index Deck - 2024-06-11 - 11.53.54.HSD5000
B4 — B4_6kfull
D4 — E4_6kfull
H4 — C5_6kfull
B5 — E5_6k full
H5 — E6_6k full
B6 — H6_6k full
4 [File 2] LP081 Run2 Index Bench Ctrls - 2024-06-10 - 16.19.51.HSD5000
D1 — 3_Ctrl_6kfull
E1 — 4_Ctrl_6kfull

20K full



4 [File 1] LP081 Run2 Index Deck - 2024-06-11 - 11.53.54.HSD5000
C4 — D4_20kfull
F4 — H4_20kfull
A5 — D5_20k full
C5 — F5_20k full
G5 — D6_20k full
A6 — G6_20k full
4 [File 2] LP081 Run2 Index Bench Ctrls - 2024-06-10 - 16.19.51.HSD5000
F1 — 5_Ctrl_20kfull
G1 — 6_Ctrl_20kfull

Figure legend: At the end of Index PCR purification, libraries were quantified via Qubit and normalized to 1 ng/μL prior to running on TapeStation for size analysis. All cell inputs prepped on the Hamilton™ Microlab™ NGS STAR ARW were comparable in bp size distribution against manual bench controls.

WTA NGS data from single cells in PBMC samples

Sequencing QC metrics




Figure legend: Three runs of 18 libraries (six each of 6K subsampled beads, 6K full sampled beads and 20K full sampled beads from the same pool of cellular transcripts captured on the BD Rhapsody™ Enhanced Cell Capture Beads) were prepared on the Hamilton™ Microlab™ NGS STAR ARW along with five manual bench controls of each condition. After sequencing, all automated and manual libraries were down sampled to between 7,600–8,000 raw reads per cell. **A)** The percentage of reads identified with a cell label UMI and align uniquely generated from libraries prepped on the Hamilton™ Microlab™ NGS STAR ARW are comparable to manually prepped libraries. **B)** Percentage of mapped reads that align to a known gene generated on the automated system are like those libraries generated manually.

WTA NGS data from single cells in PBMC samples

Cells retrieved by sequencing




Figure legend: Putative cell detection is similar between automated library preparation and manual library preparation. The expected cell inputs of 6,000 cells and 20,000 cells were approximately detected.

WTA NGS data from single cells in PBMC samples

Assay sensitivity




Figure legend: Three runs of 18 libraries (six each of 6K subsampled beads, 6K full sampled beads and 20K full sampled beads from the same pool of cellular transcripts captured on the BD Rhapsody™ Enhanced Cell Capture Beads) were prepared on the Hamilton™ Microlab™ NGS STAR ARW along with five manual bench controls of each condition. Shown are sequenced libraries of 15 libraries total for 6K sub and 6K full each and nine libraries total for 20K full sampled beads. After sequencing, all automated and manual libraries were down sampled to between 7,600–8,000 raw reads per cell. **A)** Median genes per cell generated from libraries prepped on the Hamilton™ Microlab™ NGS STAR ARW are comparable to manually prepped libraries. **B)** Median UMI counts per cell generated on the automated system are like those libraries generated manually.

WTA NGS data from single cells in PBMC samples

Assay sensitivity




Figure legend: Gene expression correlation plots for 6K subsampled, 6K full sampled and 20K full sampled beads show a correlation greater than 0.99 when automated library prepped samples are compared to manually prepared samples. Here all genes containing greater than 50 total molecules are analyzed, where 6K subsampled library contains approximately 18,000 genes and both 6K and 20K full libraries contain approximately 14,000 genes when sequenced at 8,000 raw reads per cell.

WTA NGS data from single cells in PBMC samples

Cell type identification




Figure legend: Similar percentages of cell types were identified with automated libraries prepped on the Hamilton™ Microlab™ NGS STAR ARW compared to manual bench controls.

WTA NGS data from single cells in PBMC samples

Cell type identification




Figure legend: UMAP of 6K subsamples beads shows cell types were identified and no batch effects between manual and automated library prep.

WTA NGS data from single cells in PBMC samples

Cell type identification




Figure legend: UMAP of 6K full beads shows cell types were identified and no batch effects between manual and automated library prep.

WTA NGS data from single cells in PBMC samples

Cell type identification




Figure legend: UMAP of 20K full beads shows cell types were identified and no batch effects between manual and automated library prep.

WTA NGS data from single cells in PBMC samples

Reproducibility across instruments




Figure legend: Three runs of 18 libraries (six each of 6K subsampled beads, 6K full sampled beads and 20K full sampled beads from the same pool of cellular transcripts captured on the BD Rhapsody™ Enhanced Cell Capture Beads) were prepared on the Hamilton™ Microlab™ NGS STAR ARW 1 and 2. Shown are sequenced libraries of 15 libraries total for 6K sub and 6K full each and nine libraries total for 20K full sampled beads. After sequencing, all automated and manual libraries were down sampled to between 7,600–8,000 raw reads per cell. **A)** Median genes per cell generated from libraries prepped on the Hamilton™ Microlab™ NGS STAR ARW 1 are comparable to Hamilton™ Microlab™ NGS STAR ARW 2. **B)** Median UMI counts per cell generated on the automated system are comparable across different instruments.

WTA NGS data from single cells in PBMC samples

Reproducibility across instruments




Figure legend: Wells across instruments with the same cell input were compared on UMAP and no batch effects were observed between Hamilton™ Microlab™ NGS STAR ARW 1 and Hamilton™ Microlab™ NGS STAR ARW 2.

WTA NGS data from single cells in PBMC samples

Reproducibility across instruments




Figure legend: Wells across instruments with the same cell input were compared on a gene expression correlation and the R^2 is greater than 0.99, indicating there was no difference between Hamilton™ Microlab™ NGS STAR ARW 1 and Hamilton™ Microlab™ NGS STAR ARW 2.

100,000 cell input

100,000 cell input experimental design for automated library prep:




Figure legend: 100,000 PBMCs each were loaded into Lanes A-D on two BD Rhapsody™ HT Xpress Cartridge, for a total of 8 lanes. After retrieval, approximately 68,000 to 75,000 cells per lane were retrieved for library preparation. Cartridge 1 (lanes B, C, D) were set for automated library preparation for Run 1 and Cartridge 2 (lanes B, C, D) were set for automated library preparation for Run 2. Lane A (Cartridge 1) and lane A (Cartridge 2) were prepared manually on bench.

BD Rhapsody™ Scanner metrics

Sample (Cartridge-Lane)	1-A	1-B	1-C	1-D	2-A	2-B	2-C	2-D
Number of wells with viable cells and a bead	71242	71651	72884	73014	74932	68851	73459	74473
Number of viable cells captured in wells with a bead	84391	84862	86916	87097	89335	81561	86960	88723
Cell multiplet rate	16.7%	16.7%	17.4%	17.5%	17.4%	16.7%	16.7%	17.3%
Cell retention rate	99.4%	98.9%	98.3%	99.5%	97.3%	95.2%	96.4%	97.9%
Number of retrieved beads with 1+ viable cells	71028	71436	72665	72795	74707	68644	73239	74250
Run	Manual	Run 1	Run 1	Run 1	Run 2	Run 2	Manual	Run 2
Library Prep Location	Ctrl-1	E4	A4	B4	A4	B4	Ctrl-2	E4

Figure legend: 100,000 PBMCs each were loaded into Lanes A-D on two BD Rhapsody™ HT Xpress Cartridge, for a total of 8 lanes. After retrieval, approximately 68,000 to 75,000 cells per lane were retrieved for library preparation. Cartridge 1 (lanes B, C, D) were set for automated library preparation for Run 1 and Cartridge 2 (lanes A, B, D) were set for automated library preparation for Run 2. Lane A (Cartridge 1) and lane C (Cartridge 2) were prepared manually on bench.

WTA NGS data from single cells in PBMC samples

Library DNA yield and QC metrics




Figure legend: Two runs of 3 libraries each of 100K full sampled beads of cell-captured BD Rhapsody™ Enhanced Cell Capture Beads were prepared on the Hamilton™ Microlab™ NGS STAR ARW for a total of 6 automated libraries along with two manual bench controls. **A)** Final quality control of indexed libraries show automated libraries were comparable to manual bench controls for Qubit concentration (ng/μL). **B)** Final quality control by TapeStation analysis shows similar average bp size distribution between 200–1,000bp.

WTA NGS data from single cells in PBMC samples

Final library QC TapeStation traces




Figure legend: At the end of Index PCR purification, libraries were quantified via Qubit and normalized to 1 ng/ μ L prior to running on TapeStation for size analysis. All libraries prepped on the Hamilton™ Microlab™ NGS STAR ARW were comparable in bp size distribution against manual bench controls.

WTA NGS data from single cells in PBMC samples

Sequencing QC metrics




Figure legend: Two runs of 3 libraries each of 100K full sampled beads of cell-captured BD Rhapsody™ Enhanced Cell Capture Beads were prepared on the Hamilton™ Microlab™ NGS STAR ARW for a total of 6 automated libraries along with two manual bench controls. After sequencing, all automated and manual libraries were down sampled to approximately 10,500 raw reads per cell. **A)** The percentage of reads identified with a cell label UMI and align uniquely generated from libraries prepped on the Hamilton™ Microlab™ NGS STAR ARW are comparable to manually prepped libraries. **B)** Percentage of mapped reads that align to a known gene generated on the automated system are like those libraries generated manually.

WTA NGS data from single cells in PBMC samples

Cell retrieved by sequencing




Figure legend: Putative cell detection, automated library preparation vs. manual library preparation. The expected cell inputs of 100,000 cells were approximately detected.

WTA NGS data from single cells in PBMC samples

Assay sensitivity




Figure legend: Two runs of 3 libraries each of 100K full sampled beads of cell-captured BD Rhapsody™ Enhanced Cell Capture Beads were prepared on the Hamilton™ Microlab™ NGS STAR ARW for a total of 6 automated libraries along with two manual bench controls. After sequencing, all automated and manual libraries were down sampled to approximately 10,500 raw reads per cell. **A)** Median genes per cell generated from libraries prepped on the Hamilton™ Microlab™ NGS STAR ARW are comparable to manually prepped libraries. **B)** Median UMI counts per cell generated on the automated system are like those libraries generated manually.

WTA NGS data from single cells in PBMC samples

Assay sensitivity




Figure legend: Gene expression correlation plot 100K full sampled beads show a correlation greater than 0.99 when automated library prepped samples are compared to manually prepared samples. Here all genes containing greater than 50 total molecules are analyzed, each library contains approximately 23,360 genes when sequenced at 10,500 raw reads per cell.

WTA NGS data from single cells in PBMC samples

Cell type identification




Figure legend: Similar percentages of cell types were identified with automated libraries prepped on the Hamilton™ Microlab™ NGS STAR ARW compared to manual bench controls.

WTA NGS data from single cells in PBMC samples

Cell type identification




Figure legend: UMAP of 100K full beads shows cell types were identified, and no batch effects were seen between manual and automated library prep.

Hamilton™ Microlab™ NGS STAR ARW pricing

NGS STAR pricing (NA Territory)

with MOA (MPH and ODTC add-on) deck configuration



Base Items	List Price
------------	------------

NGS STAR with 8x 1 mL Channels and iSWAP (MOA Deck)	\$160,055
NGS STAR with 8x 1mL Channels, 96 MPH, and iSWAP (MOA Deck)	\$202,233

Options	List Price
---------	------------

96 MPH Upgrade	\$42,178
On-Deck Thermal Cycler for NGS STAR (MOA Deck)	\$16,022
On-Deck Thermal Cycler for NGS STAR, 96 MPH (MOA Deck)	\$16,022

Other required materials

Ordered separately

Manufacturer	Cat. No.	Description
Hamilton	194051	60-mL reagent reservoir
Hamilton	235948	50- μ L filtered tips
Hamilton	235903	300- μ L filtered tips
Hamilton	235905	1,000- μ L filtered tips
Beckman Coulter	A63880	Agencourt AMPure XP Magnetic Beads
Fisher Scientific	BP2818100	Absolute ethyl alcohol
Invitrogen	10-977-015	Nuclease-free water
Bio-Rad	HSP9601	Hard-Shell™ 96 Well PCR Plates (low profile, thin wall, skirted, white/clear)
Thermo Fisher Scientific	AB0859	Abgene™ 96 Well 0.8 mL Polypropylene DeepWell™ Storage Plate (midi plate)

Support

Support

Supporting you with your BD[®] OMICS-One XT Automation Kit experiments



Getting help from single-cell experts

Visit us at scomix.bd.com to view our resource library, learning center and FAQs or to file a ticket for help or contact your local BD application specialist



In need of technical support

BD technical service support is here to help with assay support. Contact us **by phone** at 877.232.8995, prompt 2, then prompt 2, 6:30 a.m. to 5 p.m. PT or by **email** at ResearchApplications@bd.com.



Ordering BD[®] OMICS-One XT Assay Kits

To request a quote or place an order, visit bdbiosciences.com or contact your local BD sales representative

Support

Supporting you with your Hamilton™ Microlab™ NGS STAR ARW (NA Territory)



Getting help from automation experts

Hamilton technical support is available for your automation questions. Contact **by phone** at 1.775.858.3000, **by email** at tech@hamiltoncompany.com or contact your local Hamilton application specialist



In need of service support

Hamilton service support is here to help with instrument support. Contact **by phone** at 1.800.527.5269 or **by email** at roboticservice@hamiltoncompany.com.



Ordering the Hamilton™ Microlab™ NGS STAR ARW

To request a quote or place and order, **visit** hamiltoncompany.com/robotics, **by phone** at 1.800.646.2871, **by email** at roboticsales@hamiltoncompany.com or contact your local Hamilton sales representative



BD flow cytometers are Class 1 Laser Products.

For Research Use Only. Not for use in diagnostic or therapeutic procedures.

BD, the BD Logo and BD Rhapsody are trademarks of Becton, Dickinson and Company or its affiliates. All other trademarks are the property of their respective owners. © 2024 BD. All rights reserved. BD-135535 (v2.0) 1024