BD Rhapsody™ System Mouse TCR/BCR Next and BD OMICS-One™ WTA Next

Library Preparation Protocol

23-24999(01) For Research Use Only 2025-10

Copyrights

No part of this publication may be reproduced, transmitted, transcribed, stored in retrieval systems, or translated into any language or computer language, in any form or by any means: electronic, mechanical, magnetic, optical, chemical, manual, or otherwise, without prior written permission from BD.

The information in this guide is subject to change without notice. BD reserves the right to change its products and services at any time. Although this guide has been prepared with every precaution to ensure accuracy, BD assumes no liability for any errors or omissions, nor for any damages resulting from the application or use of this information. BD welcomes customer input on corrections and suggestions for improvement.

Patents and Trademarks

For US patents that may apply, see bd.com/patents.

BD, the BD Logo, BD Rhapsody and Omics-One are trademarks of Becton, Dickinson and Company or its affiliates. All other trademarks are the property of their respective owners. © 2025 BD. All rights reserved.

Regulatory information

For Research Use Only. Not for use in diagnostic or therapeutic procedures.

History

| Revision | Date | Change made |
|--------------|---------|------------------|
| 23-24999(01) | 2025-10 | Initial release. |

Contents

| Introduction | 5 |
|---|----|
| Symbols | 5 |
| Protocol kits | |
| Workflows | 10 |
| On bead workflow | |
| WTA library amplification workflow | |
| TCR/BCR library amplification workflow | |
| Required and recommended materials | 14 |
| · | |
| Required reagents Recommended consumables | |
| Equipment | |
| Best practices | |
| Safety information | |
| Time considerations | |
| Ducasalius | 10 |
| Procedure | 19 |
| 1. Reverse transcription, template switching, and Exonuclease I treatment | 20 |
| 1.1 cDNA Synthesis and Template Switching | 20 |
| 1.2 Denaturation and Self-hybridization | |
| 1.3 TCR/BCR Extension | |
| 1.4 Exonuclease I Treatment | 27 |
| 2. WTA library amplification | 29 |
| 2.1 WTA Random Priming and Extension (RPE) | 29 |
| 2.2 WTA RPE cleanup | 33 |
| 2.3 WTA RPE PCR | |
| 2.4 WTA RPE PCR cleanup and quantification | |
| 2.5 WTA index PCR | |
| 2.6 WTA index PCR cleanup and quality check | 44 |
| 3. TCR/BCR library amplification | 47 |
| 3.1 TCR/BCR PCR1 | 47 |
| 3.2 TCR/BCR PCR1 cleanup | |
| 3.3 TCR/BCR PCR2 | |
| 3.4 TCR/BCR PCR2 cleanup | 55 |
| 3.5 TCR/BCR RPE | |
| 3.6 TCR/BCR RPE cleanup | |
| 3.7 TCR/BCR index PCR 3.8 TCR/BCR index PCR cleanup and quality check | |
| · · · | |
| Sequencing | 67 |
| Read requirements for libraries | |
| Pooling libraries for sequencing | |
| Sequencing recommendations | |
| Seguencing analysis nineline | 70 |

| Troubleshooting | |
|--|----|
| Appendix | 73 |
| Oligonucleotides in BD OMICS-One™ WTA Next Amplification Kit | 73 |
| Mouse T cell PCR1 primers | |
| Mouse T cell PCR2 primers | 77 |
| Mouse B cell PCR1 primers | 78 |
| Mouse B cell PCR2 primers | 78 |
| Contact Information | 79 |

Introduction

This protocol enables high throughput single-cell transcriptome alongside TCR and BCR profiling of individual cells captured on the BD Rhapsody™ system, providing instructions for amplifying Illumina-compatible singlecell barcoded mRNA, TCR, and BCR libraries.

cDNA is encoded on BD Rhapsody™ Enhanced Cell Capture Beads using both the 3' and 5' ends of transcripts as templates. Whole transcriptome library is generated directly from the beads using a random priming approach, followed by an index polymerase chain reaction (PCR) step. TCR and BCR libraries are amplified from cDNA on bead using a two-step nested PCR, followed by additional random priming to capture complementarity determining regions (CDR) 1–3 and framework regions (FR) 1–4.

Symbols

The following symbols are used in this guide:

| Symbol | Description |
|----------|---|
| <u>^</u> | Important information for maintaining measurement accuracy or data integrity. |
| | Noteworthy information. |
| STOP | Procedural stopping point. |

Protocol kits

Before you begin, ensure that you have the correct kits for this protocol. Matching cap colors indicate you have the correct kit, along with the catalog numbers found in the Required and recommended materials (page 14) section.

| | BD Rhαpsody™ cDNA Kit | |
|------------|--------------------------|----------|
| Cap Color | Name | Quantity |
| | RT buffer | 1 |
| | RT 0.1M DTT | 1 |
| | Reverse transcriptase | 1 |
| | dNTP | 1 |
| | RNase inhibitor | 1 |
| • | Bead RT/PCR enhancer | 1 |
| | 10X Exonuclease I buffer | 1 |
| | Exonuclease I | 1 |
| \bigcirc | Nuclease-free water | 2 |
| | Bead Resuspension buffer | 1 |

| Cap Color | Name | Quantity |
|------------|--|----------|
| | BD Rhapsody™ HT Enhanced Cell Capture Beads v3 | 4 |
| \bigcirc | Sample buffer | 1 |
| \bigcirc | Cartridge wash buffer 1 | 1 |
| \bigcirc | Cartridge wash buffer 2 | 1 |
| \bigcirc | Lysis buffer | 4 |
| \bigcirc | Bead wash buffer | 1 |
| \bigcirc | Waste collection container | |
| \bigcirc | 1M DTT | 1 |

| | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 |
|---|---|---|---|---|---|---|---|---|----|------|
| Α | | | | | | | | | | |
| В | | | | | | | | | | |
| С | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | | |
| D | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | | |
| Ε | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | Em | ipty |

| BD OMICS-One™ WTA Next Amplificαtion Kit | | | | | |
|--|--|---------------------|------------------|--|--|
| Cap Color | Name | Part Number | Vial Placement | | |
| | BD OMICS-One™ Nuclease-Free Water | 51-9025552 | A1–A4 | | |
| | BD OMICS-One™ WTA Extension Buffer | 51-9025488 | A5 | | |
| | BD OMICS-One™ WTA Extension Primer | 51-9025467 | A6 | | |
| | BD OMICS-One™ dNTP Mixture | 51-9025491 | A7 | | |
| | BD OMICS-One™ Bead RT/PCR Enhancer | 51-9025495 | A8 | | |
| | BD OMICS-One™ WTA Extension Enzyme | 51-9025499 | A9 | | |
| | BD OMICS-One™ AbSeq Primer | 51-9025468 | A10 | | |
| | BD OMICS-One™ PCR Master Mix | 51-9025466 | B1 | | |
| | BD OMICS-One™ Universal Oligo | 51-9025553 | В2 | | |
| | BD OMICS-One™ WTA Amplification Primer | 51-9025469 | В3 | | |
| | BD OMICS-One™ Elution Buffer | 51-9025554 | B4-B8 | | |
| | BD OMICS-One™ Sample Tag PCR1 Primer | 51-9025470 | В9 | | |
| | BD OMICS-One™ Sample Tag PCR2 Primer | 51-9025471 | B10 | | |
| | BD OMICS-One™ Bead Resuspension Buffer | 51-9025555 | C9, C10, D9, D10 | | |
| | BD OMICS-One™ Library Forward Primer 1–8 | See Part numbers | C1–C8 | | |
| | BD OMICS-One™ WTA Library Reverse Primer 1–8 | for primers in rows | D1-D8 | | |
| \bigcirc | BD OMICS-One™ Multiomic Library Reverse Primer 1–8 | C–E (page 8) | E1–E8 | | |


Part numbers for primers in rows C–E

| Name | Part Number |
|--|-------------|
| BD OMICS-One™ Library Forward Primer 1 | 51-9025472 |
| BD OMICS-One™ Library Forward Primer 2 | 51-9025473 |
| BD OMICS-One™ Library Forward Primer 3 | 51-9025474 |
| BD OMICS-One™ Library Forward Primer 4 | 51-9025475 |
| BD OMICS-One™ Library Forward Primer 5 | 51-9025476 |
| BD OMICS-One™ Library Forward Primer 6 | 51-9025477 |
| BD OMICS-One™ Library Forward Primer 7 | 51-9025478 |
| BD OMICS-One™ Library Forward Primer 8 | 51-9025479 |
| BD OMICS-One™ WTA Library Reverse Primer 1 | 51-9025480 |
| BD OMICS-One™ WTA Library Reverse Primer 2 | 51-9025600 |
| BD OMICS-One™ WTA Library Reverse Primer 3 | 51-9025482 |
| BD OMICS-One™ WTA Library Reverse Primer 4 | 51-9025483 |
| BD OMICS-One™ WTA Library Reverse Primer 5 | 51-9025484 |
| BD OMICS-One™ WTA Library Reverse Primer 6 | 51-9025485 |
| BD OMICS-One™ WTA Library Reverse Primer 7 | 51-9025486 |
| BD OMICS-One™ WTA Library Reverse Primer 8 | 51-9025487 |
| BD OMICS-One™ Multiomic Library Reverse Primer 1 | 51-9025489 |
| BD OMICS-One™ Multiomic Library Reverse Primer 2 | 51-9025490 |
| BD OMICS-One™ Multiomic Library Reverse Primer 3 | 51-9025492 |
| BD OMICS-One™ Multiomic Library Reverse Primer 4 | 51-9025493 |
| BD OMICS-One™ Multiomic Library Reverse Primer 5 | 51-9025494 |
| BD OMICS-One™ Multiomic Library Reverse Primer 6 | 51-9025496 |
| BD OMICS-One™ Multiomic Library Reverse Primer 7 | 51-9025497 |
| BD OMICS-One™ Multiomic Library Reverse Primer 8 | 51-9025498 |

| | BD Rhαpsody™ Mouse TCR/BCR Next Amp | mileution Rit |
|------------|-------------------------------------|---------------|
| Cap Color | Name | Quantity |
| \bigcirc | TCR/BCR extension primers | 1 |
| \bigcirc | TCR/BCR extension buffer | 1 |
| \bigcirc | TCR/BCR extension enzyme | 1 |
| | 10 mM dNTP | 2 |
| | Nuclease-free water | 2 |
| • | Bead RT/PCR enhancer | 1 |
| | TSO Next | 1 |
| | TCR N1 primer - Mouse | 1 |
| | TCR N2 primer - Mouse | 1 |
| | BCR N1 primer - Mouse | 1 |
| | BCR N2 primer - Mouse | 1 |
| \bigcirc | PCR master mix | 1 |
| | TCR/BCR universal oligo N1 | 1 |
| | TCR/BCR universal oligo N2 | 1 |
| | Elution buffer | 2 |
| | 1M MgCl ₂ | 1 |
| | Hybridization buffer | 4 |

Workflows

On bead workflow




Denaturation Supernatant:

The mRNA template is denatured off the bead. Discard the supernatant.



Self-hybridization

Resuspend the beads in pre-warmed hybridization buffer, then gradually cool down to allow the cDNA self-hybridized on bead by CC.




VDJ region is flipped and captured by bead CC strand

Extension

Copy cell label, UMI, TCR/BCR universal oligo

Exo I


Remove unused oligo capture sequences



Final bead layout

3' sequencing template (blue): mRNA library


5' sequencing template (green): TCR/BCR library



WTA library amplification workflow

WTA RPE

Random priming on bead



WTA extension primer :

GGCTCGGAGATGTGTATAAGAGACAGNNNNNNNNN

Denature off the RPE product



WTA RPE PCR Amplify RPE product




WTA Index PCR Add Illumina adapters and indices



TCR/BCR library amplification workflow

PCR 1 Step 1:

TCR/BCR universal primer and Mouse TCR/BCR N1 primer copy target region from bead



PCR 2:

TCR/BCR universal N2 primer adds sequencing handle; Mouse TCR N2 or BCR N2 primer for nested PCR enrichment



RPE for fragmentation:

Random priming is used for fragmentation for sequencing, full-length VDJ is assembled bioinformatically. TCR and BCR RPE are done separately



Index PCR:

Add Illumina adapters and indices. TCR and BCR Index are done separately



Required and recommended materials

Required reagents

Store the reagents at the storage temperature specified on the label.

| Material | Supplier | Catalog no. |
|--|-----------------|-------------|
| BD Rhapsody™ cDNA Kit ^a | BD Biosciences | 633773 |
| BD Rhapsody™ Enhanced Cartridge Reagent Kit V3 ^{a,b} | BD Biosciences | 667052 |
| BD OMICS-One™ WTA Next Amplification Kit | BD Biosciences | 572620 |
| BD Rhapsody™ Mouse TCR/BCR Next Amplification Kit ^a | BD Biosciences | 667059 |
| AMPure® XP magnetic beads | Beckman Coulter | A63880 |
| 100% ethyl alcohol, molecular biology grade | Major supplier | _ |
| Nuclease-free water | Major supplier | _ |

a. For processing more than four libraries, two orders of this catalog number are required.

b. The Enhanced Cartridge Reagent Kit V3 must be used before beginning this protocol.

Recommended consumables

| Material | Supplier | Part number/Catalog no. |
|---|--------------------------|-------------------------|
| Pipettes (P10, P20, P200, P1000) | Major supplier | - |
| Low-retention, filtered pipette tips | Major supplier | - |
| 0.2-mL PCR 8-strip tubes | Major supplier | - |
| Axygen [®] 96–Well PCR Microplates ^a Or, | Corning | PCR96HSC |
| MicroAmp™ Optical 96–Reaction Plate ^α | Thermo Fisher Scientific | N8010560 |
| MicroAmp™ Clear Adhesive Film ^a | Thermo Fisher Scientific | 4306311 |
| 15-mL conical tube | Major supplier | - |
| DNA LoBind [®] Tubes, 1.5 mL | Eppendorf | 0030108051 |
| DNA LoBind [®] Tubes, 5.0 mL | Eppendorf | 0030108310 |
| Qubit™ Assay Tubes | Thermo Fisher Scientific | Q32856 |
| Qubit™ dsDNA HS Assay Kit | Thermo Fisher Scientific | Q32851 |
| Agilent High Sensitivity DNA Kit Or, | Agilent | 5067-4626 |
| Agilent High Sensitivity D1000 ScreenTape Agilent High Sensitivity D1000 Reagents Or, | Agilent Agilent | 5067-5584 5067-5585 |
| Agilent High Sensitivity D5000 ScreenTape Agilent High Sensitivity D5000 Reagents | Agilent Agilent | 5067-5592 5067-5593 |

 $[\]hbox{a. Recommended for processing high throughput library preparation workflows.}\\$

Equipment

| Material | Supplier | Catalog no. |
|---|--------------------------|-------------|
| Microcentrifuge for 1.5–2.0-mL tubes | Major supplier | _ |
| Microcentrifuge for 0.2-mL tubes | Major supplier | _ |
| Vortexer | Major supplier | _ |
| Digital timer | Major supplier | _ |
| Eppendorf ThermoMixer [®] C | Eppendorf | 5382000023 |
| 6-tube magnetic separation rack for 1.5-mL tubes | New England Biolabs | S1506S |
| Or, 12-tube magnetic separation rack ^a Or, | New England Biolabs | S1509S |
| Invitrogen™ DynaMag™-2 magnet ^a | Thermo Fisher Scientific | 12321D |
| Low-profile magnetic separation stand for 0.2 mL, 8-strip tubes | V&P Scientific, Inc. | VP772F4-1 |
| Magnetic Stand-96 ^b | Thermo Fisher Scientific | AM10027 |
| Qubit™ 3.0 Fluorometer or similar | Thermo Fisher Scientific | Q33216 |
| Agilent [®] 2100 Bioanalyzer | Agilent Technologies | G2940CA |
| Or, Agilent [®] 4200 TapeStation System | Agilent Technologies | G2991AA |
| Heat block | Major supplier | _ |

a. Recommended for processing greater than six samples.

Best practices



The BD Rhapsody™ Enhanced Cartridge Reagent Kit V3 (Catalog no. 667052) must be used for this protocol. The BD Rhapsody™ Mouse TCR/BCR Next Amplification Kit (Catalog no. 667059) is not compatible with the BD Rhapsody™ Enhanced Cartridge Reagent Kit (Catalog no. 664887).

Cell capture

- Ensure that the intended total cell load is 7,500–20,000. Cell loads outside this recommended range might require protocol optimization and might yield suboptimal results.
- For best results, ensure that cells have high viability before proceeding with cell capture.

b. Recommended for processing high throughput library preparation workflows.

Bead handling

 When working with BD Rhapsody™ Enhanced Cell Capture Beads, use low-retention filtered tips and LoBind® Tubes.



Never vortex the beads. Pipet-mix only.

Store BD Rhapsody™ Enhanced Cell Capture Beads at 4 °C.



Do not freeze.

• Bring Agencourt AMPure® XP magnetic beads to room temperature (15–25 °C) before use. See the AMPure® XP User's Guide for information.

Master mix preparation

- Thaw reagents (except for enzymes) at room temperature.
- Keep enzymes at -25 °C to -15 °C until ready for use.
- Return reagents to correct storage temperature as soon as possible after preparing the master mix.

Denaturation and self-hybridization

- Remove supernatant promptly after 95 °C denaturation step (≤30 seconds after placing on magnet).
- Ensure that hybridization buffer is preheated at 80 °C for at least 20 minutes before resuspending beads in step 9 of 1.2 Denaturation and Self-hybridization (page 23).



Using cold or room temperature hybridization buffer might negatively impact self-hybridization efficiency.

Supernatant handling

- Read the protocol carefully before beginning each section. Note which steps require you to keep supernatant to avoid accidentally discarding required products.
- Remove supernatants without disturbing AMPure® XP beads.
- Make and use fresh 80% ethyl alcohol within 24 hours. Adjust the volume of 80% ethyl alcohol depending on the number of libraries.


Bead amplification

- Do not proceed to thermal cycling until each tube is gently mixed by pipette to ensure uniform bead suspension. Start the thermocycler program immediately after mixing.
- Save beads after the first amplification step (WTA Random Priming and Extension (RPE) (page 29)). They must be used again for the second bead amplification step (TCR/BCR PCR1 (page 47)).

Safety information

For safety information, refer to the BD Rhapsody™ HT Single-Cell Analysis System Extended-Lysis Single-Cell Capture and cDNA Synthesis Protocol (doc ID: 23-24984) or BD Rhapsody™ HT Xpress System Extended-Lysis Single-Cell Capture and cDNA Synthesis Protocol (doc ID: 23-24983).

Time considerations



Procedure

Perform the experiment on the BD Rhapsody™ Single-Cell Analysis system following either the:

 BD Rhapsody™ HT Single-Cell Analysis Extended-Lysis System Single-Cell Capture and cDNA Synthesis Protocol (doc ID: 23-24984)



After the "Washing BD Rhapsody™ Enhanced Cell Capture Beads" section and follow this protocol from Reverse transcription, template switching, and Exonuclease I treatment (page 20) and subsequent steps.

or

BD Rhapsody™ HT Xpress System Extended-Lysis Single-Cell Capture and cDNA Synthesis Protocol (doc ID: 23-24983)



After the "Washing BD Rhapsody™ Enhanced Cell Capture Beads" section and follow this protocol from Reverse transcription, template switching, and Exonuclease I treatment (page 20) and subsequent steps.



The BD Rhapsody™ Enhanced Cartridge Reagent Kit V3 (Catalog no. 667052) must be used for this protocol.

Ensure that the intended total cell load is between 7,500-20,000 single cells for this protocol. Cell load below or above this recommended range might not be suitable for current protocol configuration. Then proceed as described in the following procedure.

1. Reverse transcription, template switching, and Exonuclease I treatment

1.1 cDNA Synthesis and Template Switching

Summary:

- Prepare cDNA mixture
- cDNA synthesis
- Add Template Switch Oligo (TSO Next)

| Item BD Part Number Preparation and Ha | | Preparation and Handling | Storage | | | | |
|--|--------------------------------------|--------------------------|--|--------|--|--|--|
| Equili | Equilibrate to room temperature: | | | | | | |
| | RT buffer | 650000067 | | | | | |
| | dNTP | 650000077 | | | | | |
| | 0.1 M DTT | 650000068 | Equilibrate to room temperature 30 minutes before setting up cDNA synthesis. Centrifuge briefly. | −20 °C | | | |
| | Nuclease-free water | 650000076 | | | | | |
| • | 1M MgCl ₂ | 91-1198 | | | | | |
| Place | Place on ice: | | | | | | |
| • | Bead RT/PCR enhancer | 91-1082 | Centrifuge briefly before adding to mix. | –20 °C | | | |
| | TSO Next | 91-1295 | | | | | |
| Leave | in freezer until ready to | use: | | | | | |
| | RNase inhibitor | 650000078 | Centrifuge briefly before adding to mix. | –20 °C | | | |
| | Reverse transcriptase | 700026321 | | | | | |
| Obtai | n: | | | | | | |
| Wash | ed enhanced cell capture l | beads | Centrifuge briefly and keep on ice until ready. | 4 °C | | | |
| Ice bu | icket | | | | | | |
| 1.5-m | 1.5-mL tube magnetic rack | | | | | | |
| 1.5-m | 1.5-mL DNA LoBind [®] tubes | | | | | | |
| Set up: | | | | | | | |
| Therm | nomixer at 42 °C | | | | | | |



This section should be performed in the pre-amplification workspace.

- 1. Set a thermomixer to 42 °C.
- 2. If performing self-hybridization on the same day, set a second thermomixer to:
 - 1,200 rpm and at 80 °C for 3 minutes.
 - 1,200 rpm and at 25 °C for 1 minute.
 - 1,200 rpm and at 25 °C infinite (optional).



The thermomixer set to 80 °C will be used as a heat block to warm the hybridization buffer, and then used with programmed cooling during 1.2 Denaturation and Self-hybridization (page 23).

3. In a new 1.5-mL LoBind $^{\circledR}$ tube, pipet the following reagents.

cDNA/template switching mix

| Сар | Component | 1 librαry (μL) | 1 library with 20% overage (μL) | 4 libraries with 20% overage (μL) | 8 libraries with 20% overage (μL) |
|-----|--------------------------|----------------|------------------------------------|-----------------------------------|-----------------------------------|
| | RT buffer | 40.0 | 48.0 | 192.0 | 384.0 |
| | dNTP | 20.0 | 24.0 | 96.0 | 192.0 |
| | RT 0.1 M DTT | 10.0 | 12.0 | 48.0 | 96.0 |
| • | Bead RT/PCR enhancer | 12.0 | 14.4 | 57.6 | 115.2 |
| | RNase inhibitor | 10.0 | 12.0 | 48.0 | 96.0 |
| | Reverse transcriptase | 10.0 | 12.0 | 48.0 | 96.0 |
| | Nuclease-free water | 98.0 | 117.6 | 470.4 | 940.8 |
| | Total | 200.0 | 240.0 | 960.0 | 1920.0 |

- 4. Gently vortex mix, briefly centrifuge, and place back on ice.
- 5. Place the tube of washed BD Rhapsody™ Enhanced Cell Capture Beads on a magnet for ≥2 minutes.
- 6. Discard the supernatant.
- 7. Remove the tube from the magnet and pipet 200 μ L of cDNA mix into the beads. Pipet-mix.



Keep the prepared cDNA mix with beads on ice until the suspension is transferred in the next step.

8. Transfer the bead suspension to a new 1.5-mL LoBind[®] tube.

9. Incubate the bead suspension on the thermomixer at 1,200 rpm and 42 °C for 30 minutes.



Shaking is critical for this incubation.

10. While the bead suspension is still incubating at 1,200 rpm and 42 °C, pipet the following reagents in a new 1.5-mL LoBind® tube.



Prepare the TSO mix approximately **within 2 minutes** before the 30 minute incubation at 42 °C is finished.



Use immediately.

TSO mix

| Сар | Component | 1 librαry (μL) | 1 library with 20% overage (μL) | 4 libraries with 20% overage (μL) | 8 libraries with 20% overage (μL) |
|-----|----------------------|----------------|------------------------------------|-----------------------------------|--------------------------------------|
| | TSO Next | 6.0 | 7.2 | 28.8 | 57.6 |
| | 1M MgCl ₂ | 2.0 | 2.4 | 9.6 | 19.2 |
| | Total | 8.0 | 9.6 | 38.4 | 76.8 |

- 11. Gently vortex mix, briefly centrifuge, and keep on ice.
- 12. Add **8 μL** of TSO mix to the reaction, gently pipet-mix, and incubate on the thermomixer for another **30 minutes** at 1,200 rpm and **42 °C**.



If you are performing self-hybridization on the same day, complete steps 3 and 4 from 1.2 Denaturation and Self-hybridization (page 23) now.

- 13. Place the bead suspension on the 1.5-mL tube magnet until the solution is clear (≤1 minute). Discard the supernatant.
- 14. Remove the tube from the magnet and pipet 200 μ L of elution buffer into the tube. Pipet-mix. Place on ice.

OPTIONAL



BD Rhapsody™ Enhanced Cell Capture Beads can be stored up to 7 days at 2–8 °C after template switching.

15. If using the BD Rhapsody™ HT Single-Cell Analysis System Instrument User Guide, view the BD Rhapsody™ Scanner image analysis to see if the analysis metrics passed.

1.2 Denaturation and Self-hybridization

Summary:

- Denature mRNA
- Add hybridization buffer to hybridize TSO onto bead

| Item | | BD Part Number | Preparation and Handling | Storage | | | | |
|---------|----------------------------------|-------------------|---|---|--|--|--|--|
| Equili | Equilibrate to room temperature: | | | | | | | |
| | Hybridization buffer | 91-1199 | Equilibrate to room temperature 30 minutes before | 20.15 | | | | |
| | Elution buffer | 91-1084 | setting up hybridization. Centrifuge briefly. | _20 °C | | | | |
| Obtain: | | | | | | | | |
| Enhar | nced cell capture beads af | ter cDNA Synthesi | s and Template Switching | | | | | |
| Ice bu | ıcket | | | | | | | |
| 1.5-m | L tube magnetic rack | | | | | | | |
| 1.5-m | L DNA LoBind [®] tubes | | | | | | | |
| Set u | Set up: | | | | | | | |
| Heat I | Heat block at 95 °C | | | | | | | |
| Heat I | Heat block at 80 °C (Optional) | | | | | | | |
| Therm | nomixer with self-hybridize | ation program | | Thermomixer with self-hybridization program | | | | |

- 1. Set a heat block to 95 °C.
- 2. Program a thermomixer with the self-hybridization program.
 - a. 1,200 rpm and at 80 °C for 3 minutes.
 - b. 1,200 rpm and at $25 \,^{\circ}\text{C}$ for 1 minute.
 - c. 1,200 rpm and at 25 °C infinite (optional).



If you performed cDNA synthesis on the same day, this is the same thermomixer from 1.1 cDNA Synthesis and Template Switching (page 20), and the thermomixer is already programmed.



Confirm "Time Mode" on the thermomixer is set to "Temperature Control" to ensure that the 25 °C temperature is reached before the 1 minute at 25 °C (step 2b) begins.

- 3. Prepare hybridization buffer for self-hybridization.
- 4. Aliquot **1.2 mL** hybridization buffer into a new 1.5-mL LoBind[®] tube and place the tube in the pre-heated 80 °C thermomixer (from step 2a) without shaking.
- 5. Keep the tube of hybridization buffer in the 80 °C thermomixer until ready to use, at least **20 minutes** before resuspending beads in step 9 of this section.
- 6. To denature, incubate the tube in the following order:
 - a. Pipet-mix to resuspend the beads.
 - b. Incubate the tube at 95 °C in a heat block for 5 minutes.
 - c. Immediately after the completion of the 95 $^{\circ}$ C incubation, slightly open the lid of the tube to release air pressure within the tube.
- 7. Immediately place the tube on the magnet for ≤30 seconds until clear.
- 8. Discard the supernatant.
- 9. Resuspend the beads in **1.0 mL** of pre-heated 80 °C hybridization buffer, and immediately place in the pre-programed thermomixer from step 2a. Start the program.



Incubation will take approximately 25 minutes.

 After the hybridization step, place tube on ice for at least 1 minute while TCR/BCR extension mix is being prepared.

1.3 TCR/BCR Extension

Summary:

- Prepare extension enzyme mix
- Extend TSO to copy cell label from bead

| Item BD Part Number | | BD Part Number | Preparation and Handling | Storage | | |
|---------------------|--------------------------------------|-----------------------|---|---------|--|--|
| Equili | Equilibrate to room temperature: | | | | | |
| \bigcirc | TCR/BCR extension buffer | 91-1206 | | | | |
| | dNTP | 650000077 | Equilibrate to room temperature 30 minutes before setting up extension. Centrifuge briefly. | –20 °C | | |
| | Nuclease-free water | 650000076 | | | | |
| Leave | in freezer until ready to | use: | | | | |
| \bigcirc | TCR/BCR extension enzyme | 91-1207 | Centrifuge briefly before adding to mix. | –20 °C | | |
| Obtai | n: | | | | | |
| Enhar | nced cell capture beads aft | er self-hybridization | Centrifuge briefly and keep on ice until ready. | 4 °C | | |
| Ice bu | ıcket | | | | | |
| 1.5-m | 1.5-mL tube magnetic rack | | | | | |
| 1.5-m | 1.5-mL DNA LoBind [®] tubes | | | | | |
| Set up | Set up: | | | | | |
| Therm | nomixer at 37 °C | | | | | |

- 1. Set a thermomixer to 37 °C.
- 2. Ensure all reagents other than the TCR/BCR extension enzyme are at room temperature.
- 3. In a new 1.5-mL LoBind $^{\circledR}$ tube, pipet the following reagents.

TCR/BCR extension mix

| Сар | Component | For 1 librαry (μL) | For 1 library with 20% overage (µL) | For 4 libraries with 20% overage (µL) | For 8 libraries with 20% overage (µL) |
|-----|--------------------------|-----------------------|--|---------------------------------------|---------------------------------------|
| | TCR/BCR extension buffer | 20 | 24 | 96 | 192 |
| | dNTP | 20 | 24 | 96 | 192 |
| | TCR/BCR extension enzyme | 10 | 12 | 48 | 96 |
| | Nuclease-free water | 150 | 180 | 720 | 1440 |
| | Total | 200 | 240 | 960 | 1920 |

- 4. Gently vortex mix, briefly centrifuge, and keep at room temperature.
- 5. Briefly spin the tube with the bead suspension.
- 6. Place the tube of BD Rhapsody™ Enhanced Cell Capture Beads on a magnet for ≤2 minutes. Discard the supernatant.
- 7. Remove the tubes from the magnet and resuspend using 200 μ L of TCR/BCR extension mix. Pipet-mix.
- 8. Incubate the bead suspension on a thermomixer at 1,200 rpm and 37 °C for 30 minutes.
- 9. Briefly spin the tube with the beads suspension and place the tube on ice.

1.4 Exonuclease I Treatment

Summary:

- Prepare exonuclease I enzyme mix
- Treat beads with exonuclease I
- Heat inactivation

| Item BD Part Number Preparatio | | Preparation and Handling | Storage | | | |
|--------------------------------|--------------------------------------|--------------------------|---|--------|--|--|
| Equili | Equilibrate to room temperature: | | | | | |
| | 10X exonuclease I buffer | 650000071 | | | | |
| \bigcirc | Nuclease-free water | 650000076 | Equilibrate to room temperature 30 minutes before setting up Exo-I treatment. Centrifuge briefly. | −20 °C | | |
| • | Bead resuspension buffer | 650000066 | | | | |
| Leave | in freezer until ready to | use: | | | | |
| | Exonuclease I | 650000072 | Centrifuge briefly before adding to mix. | –20 °C | | |
| Obtai | n: | | | | | |
| Enhar exten | nced cell capture beads af | ter TCR/BCR | Centrifuge briefly and keep on ice until ready. 4 ° | | | |
| Ice bu | ıcket | | | | | |
| 1.5-m | L tube magnetic rack | | | | | |
| 1.5-m | 1.5-mL DNA LoBind [®] tubes | | | | | |
| Set up | Set up: | | | | | |
| Therm | Thermomixer at 37 °C | | | | | |
| Heat I | olock at 80 °C | | | | | |

- 1. Set one thermomixer to 37 °C and a heat block to 80 °C.
- 2. In a new 1.5-mL LoBind[®] tube, pipet the following reagents.

Exonuclease I mix

| Сар | Component | For 1 library (μL) | For 1 library with 20% overage (µL) | For 4 libraries with 20% overage (µL) | For 8 libraries with 20% overage (µL) |
|-----|-----------------------------|-----------------------|--|---------------------------------------|---------------------------------------|
| | 10X exonuclease I buffer | 20 | 24 | 96 | 192 |
| | Exonuclease I | 10 | 12 | 48 | 96 |
| | Nuclease-free water | 170 | 204 | 816 | 1632 |
| | Total | 200 | 240 | 960 | 1920 |

- 3. Gently vortex-mix, briefly centrifuge, and keep at room temperature.
- 4. Place the tube of BD Rhapsody™ Enhanced Cell Capture Beads with TCR/BCR extension mix on a 1.5-mL tube magnet for ≤1 minute.
- 5. Discard the supernatant.
- 6. Remove the tube from the magnet and pipet 200 μ L exonuclease I mix into the tube. Pipet-mix.
- 7. Incubate the bead suspension on thermomixer at 1,200 rpm and 37 °C for 30 minutes.
- 8. Incubate the bead suspension in the heat block at 80 °C for 20 minutes.
- 9. Place the tube on ice for ~1 minute.
- 10. Briefly spin the tube with the bead suspension.
- 11. Place the tube on the magnet for ≤1 minute until clear. Discard the supernatant.
- 12. Remove the tube from the magnet and pipet **200** μ L of cold bead resuspension buffer into the tube. Pipet-mix.



Exonuclease I-treated beads can be stored in bead resuspension buffer at 4 °C for up to 1 year.

13. Proceed to library preparation.

2. WTA library amplification

2.1 WTA Random Priming and Extension (RPE)

Summary:

- Prepare random priming mix and extension enzyme mix
- Anneal random primers
- Extend random primers
- Denature RPE products

| Item | Item BD Part Number Preparation and Handling Ste | | Storage | | | |
|---------------------------------|--|----------------------|---|--------|--|--|
| Equili | Equilibrate to room temperature: | | | | | |
| | WTA extension buffer | 51-9025488 | | | | |
| | WTA extension primer | 51-9025467 | Equilibrate to room temperature 30 minutes before | | | |
| | dNTP mixture | 51-9025491 | setting up RPE. | −20 °C | | |
| \bigcirc | Nuclease-free water | 51-9025552 | Centrifuge briefly. | | | |
| | Elution buffer | 51-9025554 | | | | |
| Place | on ice: | | | | | |
| | Bead RT/PCR enhancer | 51-9025495 | Centrifuge briefly before adding to mix. | −20 °C | | |
| Leave | in freezer until ready to | use: | | | | |
| | WTA extension enzyme | 51-9025499 | Centrifuge briefly before adding to mix. | −20 °C | | |
| Obtai | n: | | | | | |
| Exonu | ıclease I-treated enhancec | l cell capture beads | Centrifuge briefly and keep on ice until ready. | 4 °C | | |
| Ice bu | ıcket | | | | | |
| 1.5-m | L tube magnetic rack | | | | | |
| 1.5-m | L DNA LoBind [®] tubes | | | | | |
| Set up: | | | | | | |
| Heat I | Heat block at 95 °C | | | | | |
| Thermomixer at 37 °C (Optional) | | | | | | |
| Therm | nomixer at 25 °C | | | | | |
| Progra | ammed thermomixer with | RPE program | | | | |

This section describes how to generate random priming products. First, random primers are hybridized to the cDNA on the BD Rhapsody™ Enhanced Cell Capture Beads, followed by extension with an enzyme.



Perform this procedure in the pre-amplification workspace. We recommend using a separate heat block for the $95\,^{\circ}\text{C}$ incubations.

1. Set a heat block to 95 °C and set two thermomixers to 37 °C and 25 °C, respectively.



If you are using one thermomixer, skip the 37 °C incubation in step 11b.

2. In a new 1.5-mL LoBind[®] tube, pipet the following reagents.

Random primer mix

| Сар | Component | 1 librαry (μL) | 1 library with 20% overage (μL) | 4 libraries with 20% overage (μL) | 8 libraries with 20% overage (μL) |
|-----|----------------------|----------------|------------------------------------|-----------------------------------|--------------------------------------|
| | WTA extension buffer | 20.0 | 24.0 | 96.0 | 192.0 |
| | WTA extension primer | 40.0 | 48.0 | 192.0 | 384.0 |
| | Nuclease-free water | 114.0 | 136.8 | 547.2 | 1,094.4 |
| | Total | 174.0 | 208.8 | 835.2 | 1670.4 |

- 3. Pipet-mix the random primer mix and keep at room temperature.
- 4. Briefly centrifuge the tube of Exonuclease I-treated BD Rhapsody™ Enhanced Cell Capture Beads and then complete one of the following actions.
 - If you are using a subsample of the beads, proceed to the next step.
 - If you are using the entire sample of beads, skip to step 6.
- 5. (Optional) To subsample the Exonuclease I-treated BD Rhapsody™ Enhanced Cell Capture Beads:
 - a. Based on the expected number of viable cells captured on the beads in the final bead resuspensionvolume, determine the volume of beads to subsample for sequencing.
 - b. Completely resuspend the beads by pipet-mixing, then pipet the calculated volume of the bead suspension into a new 1.5-mL LoBind $^{\textcircled{\$}}$ tube.



The remaining Exonuclease I-treated beads can be stored in bead resuspension buffer at $4 \,^{\circ}$ C for up to 1 year.



Subsample is only optional for WTA + TCR/BCR.

- 6. Place the tube on a magnet until the supernatant is clear (<2 minutes).
- 7. Remove and discard the supernatant.
- 8. Remove the tube from the magnet.

- 9. Add 174 μ L of random primer mix into the tube.
- 10. Pipet-mix 10 times to resuspend the beads.
- 11. Incubate the tube in the following order:
 - a. 95 °C in a heat block (no shaking) for 5 minutes.
 - b. Thermomixer at 1,200 rpm and at 37 °C for 5 minutes.



Optional: If you are using one thermomixer, skip the 37 °C incubation.

- c. Thermomixer at 1,200 rpm and at 25 °C for 5 minutes.
- 12. Briefly centrifuge the tube.
- 13. Place at room temperature until ready to use.
- 14. In a new 1.5-mL LoBind[®] tube, pipet the following reagents.

Extension enzyme mix

| Сар | Component | For 1 library (µL) | For 1 library with 20% overage (µL) | For 4 libraries with 20% overage (µL) | For 8 libraries with 20% overage (µL) |
|-----|-------------------------|-----------------------|--|---------------------------------------|---------------------------------------|
| | dNTP | 8.0 | 9.6 | 38.4 | 76.8 |
| | Bead RT/PCR enhancer | 12.0 | 14.4 | 57.6 | 115.2 |
| | WTA extension enzyme | 6.0 | 7.2 | 28.8 | 57.6 |
| | Total | 26.0 | 31.2 | 124.8 | 249.6 |

- 15. Pipet 26 µL of the extension enzyme mix into the sample tube containing the beads (for a total volume of 200 $\mu L)$ and keep on ice until ready.
- 16. Program the thermomixer.
 - a. 1,200 rpm and at 25 °C for 10 minutes.
 - b. 1,200 rpm and at 37 °C for 15 minutes.
 - c. 1,200 rpm and at 45 °C for 10 minutes.
 - d. 1,200 rpm and at 55 °C for 10 minutes.



Confirm "Time Mode" is set to "Time Control" before the program begins.

- 17. Place the tube of extension enzyme mix with BD Rhapsody™ Enhanced Cell Capture Beads in the programmed thermomixer. The program takes 45 minutes.
- 18. Remove the tube after the program is complete.
- 19. Place the tube on a magnet until the supernatant is clear (<2 minutes).

- 20. Remove and discard the supernatant.
- 21. Remove the tube from the magnet.
- 22. Pipet **200 μL** of elution buffer into the tube.
- 23. Pipet-mix 10 times until the beads are fully resuspended.
- 24. Place the tube on a magnet until the supernatant is clear (<2 minutes).
- 25. Remove and discard the supernatant.
- 26. Remove the tube from the magnet.
- 27. Pipet **80** μ L of elution buffer into the tube.
- 28. To denature the random priming products off the beads.
 - a. Pipet-mix 10 times to resuspend the beads.
 - b. Incubate the tube at 95 °C in a heat block for 5 minutes (no shaking).
 - c. Slightly open the lid of the tube to release air pressure within the tube.
 - d. Place the tube on ice for 1 minute.
 - e. Briefly centrifuge the tube.
 - f. Place the tube on a magnet until the supernatant is clear (<2 minutes).



SAVE SUPERNATANT AT THIS STEP. Do not discard!

- g. Transfer **80 \muL** of the supernatant (RPE product) to a new 1.5-mL LoBind $^{\circledR}$ tube.
- 29. Place the tube containing the RPE product on ice. The total volume of RPE product will be 80 μ L. Proceed to WTA RPE cleanup (page 33).
- 30. Pipet **200** µL of cold bead resuspension buffer to the tube with leftover beads. Gently resuspend the beads by pipet-mixing only. Do not vortex.
- 31. Store the beads on ice or at 4 °C in the pre-amplification workspace until needed.



These beads will be used for TCR/BCR library amplification (page 47). DO NOT THROW AWAY!

2.2 WTA RPE cleanup

Summary:

• RPE cleanup

| Item | | BD Part Number | Preparation and Handling | Storage | | | | | |
|--------------------------------------|-----------------------------------|----------------|----------------------------------|---------|--|--|--|--|--|
| Equilibrate to room temperature: | | | | | | | | | |
| | Elution buffer | 51-9025554 | Centrifuge briefly. | −20 °C | | | | | |
| AMPu | re [®] XP magnetic beads | | Manufactural | | | | | | |
| Qubit dsDNA HS Assay Kit | | | - Manufacturer's recommendations | | | | | | |
| Obtain: | | | | | | | | | |
| WTA RPE product | | | | 4 °C | | | | | |
| 1.5-mL tube magnetic rack | | | | | | | | | |
| 1.5-mL DNA LoBind [®] tubes | | | | | | | | | |
| Set up: | | | | | | | | | |
| Prepare fresh 80% ethyl alcohol | | | | | | | | | |

This section describes how to perform a single-sided AMPure[®] cleanup, which removes primer dimers and other small molecular weight byproducts. The final product is purified single-stranded DNA.



Perform the RPE purification in the pre-amplification workspace.

1. Make fresh 80% ethyl alcohol for use within 24 hours.



Adjust the volume depending on the number of samples. One sample requires 0.5 mL of 80% ethyl alcohol.

- 2. Bring AMPure[®] XP beads to room temperature. Vortex the AMPure[®] XP beads at high speed for **1 minute** until the beads are fully resuspended.
- 3. Pipet 128 μ L of AMPure[®] XP beads into the tube containing the 80 μ L of RPE product supernatant (1.6x). Pipet-mix at least 10 times, then briefly centrifuge.
- 4. Incubate at room temperature for 5 minutes.
- 5. Place the tube on the magnet for 3 minutes. Discard the supernatant.
- Keeping the tube on the magnet, gently add 250 μL of fresh 80% ethyl alcohol into the tube and incubate for 30 seconds. Discard the supernatant.
- 7. Repeat step 6 for a total of two ethyl alcohol washes.
- 8. Keeping the tube on the magnet, use a P20 pipette to remove and discard any residual supernatant from the tube.
- 9. Air-dry the beads at room temperature for 5 minutes or until the beads no longer look glossy.
- 10. Remove the tube from the magnet and resuspend the bead pellet in 80 μ L of elution buffer.
- 11. Incubate the sample at room temperature for **2 minutes**. Briefly centrifuge the tube to collect the contents at the bottom.
- 12. Place the tube on the magnet until the solution is clear, usually ~30 seconds.
- 13. Pipet the eluate (\sim 80 μ L) to a new PCR tube. This is the purified RPE product.
- 14. Keep on ice until ready to proceed with WTA RPE PCR (page 35).

2.3 WTA RPE PCR

Summary:

- Prepare RPE PCR mix
- Amplify using RPE PCR program

| Item | | BD Part Number | Preparation and Handling | Storage | | | | | |
|--------------------------------------|--------------------------|-------------------|---|---------|--|--|--|--|--|
| Equilibrate to room temperature: | | | | | | | | | |
| | Universal oligo | 51-9025553 | Cavilibrate to veces town partius 20 minutes before | −20 °C | | | | | |
| | WTA amplification primer | 51-9025469 | Equilibrate to room temperature 30 minutes before setting up RPE PCR. Centrifuge briefly. | | | | | | |
| Leave in freezer until ready to use: | | | | | | | | | |
| | PCR master mix | 51-9025466 | Centrifuge briefly before adding to mix. | −20 °C | | | | | |
| Obtain: | | | | | | | | | |
| RPE product | | | | | | | | | |
| Ice bucket | | | | | | | | | |
| 0.2-mL PCR tubes | | | | | | | | | |
| Set up: | | | | | | | | | |
| Thermocycler with RPE PCR program | | | | | | | | | |

This section describes how to generate more RPE product through PCR amplification, so that there are multiple copies of each random-primed molecule.



Perform this section in the pre-amplification workspace.

1. In the pre-amplification workspace, in a new 1.5-mL LoBind $^{\textcircled{\$}}$ tube, pipet the following components.

RPE PCR mix

| Сар | Component | For 1 librαry (μL) | For 1 library with 20% overage (µL) | For 4 libraries with 20% overage (µL) | For 8 libraries with 20% overage (µL) |
|-----|--------------------------|-----------------------|--|---------------------------------------|---------------------------------------|
| | PCR master mix | 30.0 | 36.0 | 144.0 | 288.0 |
| | Universal oligo | 6.0 | 7.2 | 28.8 | 57.6 |
| | WTA amplification primer | 6.0 | 7.2 | 28.8 | 57.6 |
| | Total | 42.0 | 50.4 | 201.6 | 403.2 |

- 2. Pipet-mix the RPE PCR mix.
- 3. Place on ice until ready to use.
- 4. Add $42~\mu L$ of the RPE PCR mix to the tube with the $80~\mu L$ of RPE product.
- 5. Pipet-mix 10 times to create the RPE PCR reaction mix.
- 6. Split the RPE PCR reaction mix into two 0.2-mL PCR tubes with $61~\mu L$ mix per tube.
- 7. Transfer any residual mix to one of the tubes.



Bring the tubes to the post-amplification workspace.

8. Run the following PCR program.

RPE PCR program

| Step | Cycles | Temperature | Time |
|-----------------|--|-------------|------------|
| Hot start | 1 | 98 °C | 45 seconds |
| Denaturation | | 98 ℃ | 15 seconds |
| Annealing | Recommended PCR cycles for resting peripheral blood mononuclear cells (PBMCs)* | 60 °C | 30 seconds |
| Extension | 7,500 – 20,000 cells: 10 cycles | 72 °C | 1 minute |
| Final extension | 1 | 72 ℃ | 2 minutes |
| Hold | 1 | 4 °C | ∞ |

^{*}Recommended number of PCR cycles might require optimization for different cell types.



The PCR can run overnight.

9. When the RPE PCR program is complete, briefly centrifuge the tubes.

2.4 WTA RPE PCR cleanup and quantification

Summary:

- RPE PCR cleanup (2 rounds)
- Quantify using Qubit Fluorometer

| Item | | BD Part Number | Preparation and Handling | Storage | |
|--------|--------------------------------------|----------------|--------------------------------|---------|--|
| Equili | brate to room temperatu | ire: | | | |
| | Elution buffer | 51-9025554 | Centrifuge briefly. | –20 °C | |
| AMPu | re [®] XP magnetic beads | | Manufacturer's recommendations | | |
| Qubit | dsDNA HS Assay Kit | | | | |
| Obtai | n: | | | | |
| RPE P | RPE PCR product | | | | |
| 1.5-m | 1.5-mL DNA LoBind [®] tubes | | | | |
| 0.2-m | 0.2-mL PCR tubes | | | | |
| 1.5-m | 1.5-mL tube magnetic rack | | | | |
| Set u | Set up: | | | | |
| Prepa | re fresh 80% ethyl alcohol | | | | |



Perform the purification in the post-amplification workspace.

- 1. Bring AMPure[®] XP beads to room temperature.
- 2. Make fresh 80% ethyl alcohol and use within 24 hours.



Adjust the volume of 80% ethyl alcohol depending on the number of samples. One sample requires 1 mL 80% ethyl alcohol.

- 3. Vortex the AMPure® XP beads until the beads are fully resuspended.
- 4. Briefly centrifuge the tubes with the RPE PCR product.
- Combine the **two** tubes of **61 \muL** RPE PCR into a new 1.5-mL LoBind[®] tube.
- Pipet-mix 10 times.
- 7. Transfer exactly **110 \muL** RPE PCR product to a new 1.5-mL LoBind[®] tube.
- Pipet 88 μ L of AMPure (0.8x) into the tube.
- 9. Pipet-mix 10 times.
- 10. Briefly centrifuge the tube.



Avoid getting AMPure® XP beads on the lid of the tube. Residual AMPure® XP beads and PCR mix buffer can negatively impact downstream results.

- 11. Incubate at room temperature for **5 minutes**.
- Place the tube on a magnet until the supernatant is clear (<5 minutes).
- Remove and discard the supernatant.
- Keeping the tube on the magnet, gently pipet 200 µL of fresh 80% ethyl alcohol into the tube.
- 15. Incubate for 30 seconds.
- 16. Remove and discard the supernatant without disturbing the beads.
- Repeat steps 14–17 once for a total of two ethyl alcohol washes.
- 18. Keeping the tube on the magnet, use a P20 pipette to remove and discard any residual supernatant from the tube.
- 19. Air-dry the beads at room temperature until the beads no longer look glossy (~3 minutes)



Do not overdry the AMPure® XP beads after the ethyl alcohol washes. Overdried beads appear cracked.

- 20. Remove the tube from the magnet.
- 21. Pipet 40 μ L of elution buffer into the tube.
- 22. Pipet-mix 10 times until the beads are fully resuspended.

- 23. Incubate at room temperature for 2 minutes.
- 24. Briefly centrifuge the tube.
- 25. Place the tube on a magnet until the supernatant is clear (~30 seconds).
- 26. Pipet the eluate (40 μ L) into a new 0.2-mL PCR strip tube.
- 27. Add **60 \muL** of water to the eluate for a final volume of **100 \muL**.



The volume must be exactly $100 \mu L$. Adjust the water volume if needed.

- 28. Pipet **80** μ L of AMPure[®] XP beads (**0.8x**) into the tube.
- 29. Pipet-mix 10 times.
- 30. Briefly centrifuge the tube.
- 31. Incubate at room temperature for 5 minutes.
- 32. Place the tube on a magnet until the supernatant is clear (<5 minutes).
- 33. Remove and discard the supernatant.
- 34. Keeping the tube on the magnet, gently pipet 200 µL of fresh 80% ethyl alcohol into the tube.
- 35. Incubate for 30 seconds.
- 36. Remove and discard the supernatant without disturbing the beads.
- 37. Repeat steps 34–36 for a total of two ethyl alcohol washes.
- 38. Keeping the tube on the magnet, use a P20 pipette to remove and discard any residual supernatant from the tube.
- 39. Air-dry the beads at room temperature until the beads no longer look glossy (~3 minutes).
- 40. Remove the tube from the magnet.
- 41. Pipet 30 μ L of elution buffer into the tube.
- 42. Pipet-mix 10 times until the beads are fully resuspended.
- 43. Incubate at room temperature for 2 minutes.
- 44. Briefly centrifuge the tube.
- 45. Place the tube on a magnet until the supernatant is clear (~30 seconds).
- 46. Pipet the eluate (30 μ L) into a new PCR strip tube.

The purified RPE PCR product is ready for Section 2.5: WTA index PCR (page 41).

47. Quantify the RPE PCR products with a Qubit™ Fluorometer using the Qubit™ dsDNA HS Assay.



Purified PCR product can be stored at -20 °C for up to 6 months.

2.5 WTA index PCR

Summary:

- Prepare WTA index PCR mix
- Amplify using WTA index PCR program

| Item | | BD Part Number | Preparation and Handling | Storage | | |
|---------|---|----------------|--|---------|--|--|
| Equili | brate to room temperature: | | | | | |
| | Forward primer 1–8 | Various | Equilibrate to room temperature 30 minutes | | | |
| | WTA reverse primer 1–8 | Various | before setting up WTA index PCR. | −20 °C | | |
| | Nuclease-free water | 51-9025552 | Centrifuge briefly. Keep on ice until ready. | | | |
| Leave | Leave in freezer until ready to use: | | | | | |
| | PCR master mix | 51-9025466 | Centrifuge briefly before adding to mix. | −20 °C | | |
| Obtai | n: | | | | | |
| Purifie | Purified RPE PCR product 4 °C | | | | | |
| Ice bu | ıcket | | | | | |
| 1.5-m | 1.5-mL DNA LoBind [®] tubes | | | | | |
| 0.2-m | 0.2-mL PCR tubes | | | | | |
| Set up | Set up: | | | | | |
| Therm | Thermocycler with WTA index PCR program | | | | | |

This section describes how to generate mRNA libraries compatible with various sequencing platforms, by adding full-length sequencing adapters and indices through PCR. We provide reagents for unique dual-indexing, with different library forward primers and reverse primers for up to 8 samples.

The same indices can be used for all library types for each lane (WTA, TCR, and BCR, for example). The libraries will be demultiplexed using the BD Rhapsody™ Sequence Analysis Pipeline. If you prefer to index each library separately, you can use combinatorial dual indexing for more index combinations.



Consult sequencing platform guidelines for low-plex pooling, to ensure the indices chosen meet the color balancing guidelines for the sequencing instrument that will be used.

1. In a new 1.5-mL LoBind $^{\textcircled{8}}$ tube, pipet the following components.

WTA index PCR mix

| Сар | Component | For 1 library (µL) | For 1 library with 20% overage (µL) | For 4 libraries with 20% overage (µL) | For 8 libraries with 20% overage (µL) |
|------------|------------------------|-----------------------|-------------------------------------|---------------------------------------|---------------------------------------|
| | PCR master mix | 12.5 | 15.0 | 60.0 | 120.0 |
| | Forward primer 1–8 | 2.5 | 3.0 | N/A | N/A |
| | WTA reverse primer 1–8 | 2.5 | 3.0 | N/A | N/A |
| \bigcirc | Nuclease-free water | 22.5 | 27.0 | 108.0 | 216.0 |
| | Total | 40.0 | 48.0 | 168.0 | 336.0 |

- 2. Pipet-mix the WTA index PCR mix.
- 3. Pipet 35 μ L into separate 0.2-mL PCR tubes for each sample.
- 4. Add 2.5 μ L of forward primer and 2.5 μ L of reverse primer to each sample.
- 5. Place on ice until ready to use.
- 6. Dilute an aliquot of the purified RPE PCR product from step 46 of WTA RPE PCR cleanup and quantification (page 38) with water to **0.5 ng/µL**.



If RPE PCR product concentration is <0.5 $\,$ ng/ μ L, adjust the number of index PCR cycles as outlined in the table.

- 7. Add 10 μ L of RPE PCR product to 40 μ L index PCR mix.
- 8. Pipet-mix 10 times.

9. Run the following PCR program.

WTA index PCR program

| Step | Cycles | Temperature | Time |
|-----------------|--|-------------|------------|
| Hot start | 1 | 98 °C | 45 seconds |
| Denaturation | RPE PCR concentration* | 98 ℃ | 15 seconds |
| Annealing | < 0.2 ng/μL: 11 cycles 0.2 ng/μL: 10 cycles | 60 ℃ | 30 seconds |
| Extension | 0.5 ng/μL: 8 cycles | 72 °C | 1 minute |
| Final extension | 1 | 72 °C | 2 minutes |
| Hold | 1 | 4 °C | ∞ |

^{*}Recommended number of PCR cycles might require optimization for different cell types.



The PCR can run overnight.

10. When the WTA index PCR program is complete, briefly centrifuge the tubes.

2.6 WTA index PCR cleanup and quality check

Summary:

- WTA index PCR cleanup
- Quality check using Qubit Fluorometer and BioAnalyzer/TapeStation

| Item | | BD Part Number | Preparation and Handling | Storage | |
|------------|---|----------------|--------------------------------|---------|--|
| Equili | Equilibrate to room temperature: | | | | |
| | Elution buffer | 51-9025554 | | 20.86 | |
| \bigcirc | Nuclease-free water | 51-9025552 | Centrifuge briefly. | –20 °C | |
| AMPu | re [®] XP magnetic beads | | | | |
| Qubit | dsDNA HS Assay Kit | | | | |
| OR | Agilent BioAnalyzer High Sensitivity Kit OR Agilent TapeStation ScreenTape and Reagents | | Manufacturer's recommendations | | |
| Obtai | n: | | | | |
| WTA | index PCR product | | | 4 °C | |
| 1.5-m | L DNA LoBind [®] tubes | | | | |
| 0.2-m | 0.2-mL PCR tubes | | | | |
| 0.2-m | 0.2-mL PCR tube magnetic rack | | | | |
| Set u | Set up: | | | | |
| Prepa | re fresh 80% ethyl alcoho | ol | | · | |

This section describes how to perform a single-sided AMPure® XP beads cleanup for sequencing. The final product is purified double-stranded DNA with full-length adapter sequences.



Perform the purification in the post-amplification workspace.

- 1. Bring AMPure[®] XP beads to room temperature.
- 2. Make fresh 80% ethyl alcohol and use within 24 hours. Adjust the volume of 80% ethyl alcohol depending on the number of samples. One sample requires 0.5 mL 80% ethyl alcohol.
- 3. Vortex the AMPure® XP beads until the beads are fully resuspended.
- 4. Add **60 \muL** of water to **50 \muL** of the WTA index PCR product.
- 5. Transfer 100 μ L of WTA index PCR product into a new 0.2-mL PCR tube.



The volume must be exactly 100 μ L.

- 6. Pipet **65** μ L of AMPure[®] XP beads (**0.65**x) into the tube.
- 7. Pipet-mix 10 times.
- 8. Briefly centrifuge the tube.
- 9. Incubate at room temperature for 5 minutes.
- 10. Place the tube on a magnet until the supernatant is clear (<5 minutes).
- 11. Remove and discard the supernatant.
- Keeping the tube on the magnet, gently pipet 200 μ L of fresh 80% ethyl alcohol into the tube.
- 13. Incubate for 30 seconds.
- 14. Remove and discard the supernatant without disturbing the beads.
- 15. Repeat steps 12–14 once for a total of two ethyl alcohol washes.
- 16. Keeping the tube on the magnet, use a P20 pipette to remove and discard any residual supernatant from the tube.
- 17. Air-dry the beads at room temperature until the beads no longer look glossy (~2 minutes).
- 18. Remove the tube from the magnet.
- 19. Pipet 30 μ L of elution buffer into the tube.
- Pipet-mix 10 times until the beads are fully resuspended.
- 21. Incubate at room temperature for 2 minutes.
- 22. Briefly centrifuge the tube.
- 23. Place the tube on the magnet until the solution is clear (~30 seconds).
- 24. Pipet the eluate (30 μ L) into a new 1.5-mL LoBind[®] tube. The purified eluate is the final sequencing library.




Purified PCR product can be stored at -20 °C for up to 6 months.

- 25. Quantify and perform quality control of the WTA index PCR product with a Qubit™ Fluorometer using the Qubit™ dsDNA HS Assay and one of the following systems:
 - Agilent 2100 BioAnalyzer using the Agilent High Sensitivity DNA Kit
 - Agilent 4200 TapeStation system using the Agilent High Sensitivity D1000 or D5000 ScreenTape assay

The expected concentration from the Qubit™ Fluorometer is >1 ng/µL.

The TapeStation trace should show a peak from \sim 200 to 2,000 bp. Refer to the representative traces in the following figures.

Figure 1 Representative TapeStation High Sensitivity D5000 trace – WTA index PCR product



3. TCR/BCR library amplification

3.1 TCR/BCR PCR1

Summary:

- Prepare TCR/BCR PCR1 mix
- Amplify using TCR/BCR PCR1 program

| Item | | BD Part Number | Preparation and Handling | Storage | | | |
|--|----------------------------------|----------------|--|---------|--|--|--|
| Equili | Equilibrate to room temperature: | | | | | | |
| 0 | PCR master mix | 91-1083 | | | | | |
| | TCR/BCR universal oligo N1 | 91-1204 | | −20 °C | | | |
| | TCR N1 primer - Mouse | 91-1212 | Equilibrate to room temperature 30 minutes before setting up TCR/BCR PCR1. Centrifuge briefly. Keep on ice | | | | |
| | BCR N1 primer - Mouse | 91-1214 | until ready. | | | | |
| • | Bead RT/PCR enhancer | 91-1082 | | | | | |
| | Nuclease-free water | 650000076 | | | | | |
| Obtai | n: | | | | | | |
| Enhanced cell capture beads after WTA RPE Denaturation | | | | | | | |
| Ice bucket | | | | | | | |
| 0.2-mL PCR tubes | | | | | | | |
| Set up | Set up: | | | | | | |
| Therm | nocycler with TCR/BCR PCF | R1 program | | | | | |

- 1. Obtain beads from step 30 of 2.1 WTA Random Priming and Extension (RPE) (page 29).
- 2. In the pre-amplification workspace, pipet the following reagents into a new 1.5-mL LoBind $^{\textcircled{8}}$ tube.

TCR/BCR PCR1 mix

| Сар | Component | For 1 library (µL) | For 1 library with 20% overage (µL) | For 4 libraries with 20% overage (µL) | For 8 libraries with 20% overage (µL) |
|------------|-------------------------------------|-----------------------|--|---------------------------------------|---------------------------------------|
| \bigcirc | PCR master mix | 100.0 | 120.0 | 480.0 | 960.0 |
| | TCR/BCR universal oligo N1 | 10.0 | 12.0 | 48.0 | 96.0 |
| • | Bead RT/PCR enhancer | 12.0 | 14.4 | 57.6 | 115.2 |
| | Mouse TCR N1 primer ^a | 2.4 | 2.9 | 11.5 | 23.0 |
| | Mouse BCR N1 primer ^a | 2.4 | 2.9 | 11.5 | 23.0 |
| | Nuclease-free water | 73.2 | 87.8 | 351.4 | 702.8 |
| | Total | 200.0 | 240.0 | 960.0 | 1920.0 |

a. If only doing TCR or BCR amplification, replace N1 primer volume with water. For example, if only doing TCR amplification, replace BCR N1 primer with water.

- 3. Gently vortex mix, briefly centrifuge, and place back on ice.
- 4. Briefly spin the tube with the bead suspension.
- 5. Place the tube of beads on a magnet for ≤1 minute.
- 6. Discard the supernatant.
- 7. Remove the tube from the magnet and resuspend the beads in **200 \muL** of TCR/BCR PCR1 mix to create the TCR/BCR PCR1 reaction mix.
- 8. Do not vortex.
- 9. Ensuring that the beads are fully resuspended, pipet **50** μ L of TCR/BCR PCR1 reaction mix with beads into each of four 0.2-mL PCR tubes.
- 10. Transfer any residual mix to one of the tubes.



Bring the TCR/BCR PCR1 reaction mix to the post-amplification workspace.

11. Run the following PCR program on the thermal cycler.

| TCR/BCR PCR1 program |
|----------------------|
|----------------------|

| Step | Cycles | Temperature | Time |
|-----------------|-------------------------------------|--------------------|------------|
| Hot start | 1 | 95 °C ^α | 3 minutes |
| Denaturation | Recommended PCR cycles ^b | 95 ℃ | 30 seconds |
| Annealing | 7,500 – 10,000 cells: 11 cycles | 60 °C | 1 minute |
| Extension | 20,000 cells: 10 cycles | 72 °C | 1 minute |
| Final extension | 1 | 72 ℃ | 5 minutes |
| Hold | 1 | 4 °C | ∞ |

a. To avoid beads settling due to prolonged incubation time on the thermal cycler before the denaturation step, it is critical to pause the instrument at 95 °C before loading the samples. Different thermal cyclers might have different pause time settings. In certain brands of thermal cyclers, however, we have observed a step-skipping error with the pause/unpause functions. To ensure that the full 3-minute denaturation is not skipped, verify that the pause/unpause functions are working correctly on your thermal cycler. To avoid the step-skipping problem, a 1-minute 95 °C pause step can be added immediately before the 3-minute 95 °C denaturation step.

12. Ramp the heated lid and heat block of the post-amplification thermal cycler to ≥95 °C by starting the thermal cycler program and then pausing it.



Do not proceed to thermal cycling until each tube is gently mixed by pipette to ensure uniform bead suspension.

13. For each 0.2-mL PCR tube, gently pipet-mix, immediately place the tube in thermal cycler, and unpause the thermal cycler program.



The PCR can run overnight, but proceed with purification within 24 hours after PCR.

- 14. After PCR, briefly centrifuge the tubes.
- 15. Put the tubes on a magnet for >30 seconds.
- 16. For each sample, collect and combine the supernatant from the four 0.2-mL PCR tubes into one new 1.5-mL LoBind[®] tube without disturbing the beads.
- 17. Discard the beads.

b. Recommended PCR cycles might need to be optimized for different cell types and cell number.

3.2 TCR/BCR PCR1 cleanup

Summary:

• TCR/BCR PCR1 cleanup

| Item | | BD Part Number | Preparation and Handling | Storage | | |
|------------|-----------------------------------|----------------|--------------------------------|---------|--|--|
| Equili | Equilibrate to room temperature: | | | | | |
| | Elution buffer | 91-1084 | Contribute by the | 20.00 | | |
| \bigcirc | Nuclease-free water | 65000076 | Centrifuge briefly. | –20 °C | | |
| AMPu | re [®] XP magnetic beads | | Manufacturer's recommendations | | | |
| Obtai | Obtain: | | | | | |
| TCR/E | SCR PCR1 product | | | 4 °C | | |
| 1.5-m | L DNA LoBind [®] tubes | | | | | |
| 0.2-m | 0.2-mL PCR tubes | | | | | |
| 1.5-m | 1.5-mL tube magnetic rack | | | | | |
| Set u | Set up: | | | | | |
| Prepa | Prepare fresh 80% ethyl alcohol | | | | | |

This section describes how to perform a single-sided AMPure® XP beads cleanup to remove primer dimers from the TCR/BCR PCR1 products. The final product is purified double-stranded DNA.



Perform the purification in the post-amplification workspace.

1. Make fresh 80% ethyl alcohol for use within 24 hours.



Adjust the volume depending on the number of samples. One sample requires 1 mL of 80% ethyl alcohol.

- 2. Bring the AMPure® XP beads to room temperature. Vortex at a high speed for 1 minute until the beads are fully resuspended.
- 3. To **200 \muL** of TCR/BCR PCR1 products, pipet **140 \muL** AMPure[®] XP bead (**0.7x**) (from step 17 in TCR/BCR PCR1 (page 47)).
- 4. Pipet-mix 10 times. Incubate at room temperature for 5 minutes.
- 5. Place the 1.5-mL LoBind[®] tube on the magnet for **5 minutes**.
- 6. Discard the supernatant.
- 7. Keeping the tube on the magnet, gently add 500 μ L of fresh 80% ethyl alcohol into the tube and incubate for 30 seconds.
- 8. Discard the supernatant.
- 9. Repeat steps 7–8 once for a total of two ethyl alcohol washes.
- 10. Keeping the tube on the magnet, use a P20 pipette to remove and discard any residual supernatant from the tube.
- 11. Air-dry the beads at room temperature for **5 minutes**.
- 12. Remove the tube from the magnet and resuspend the bead pellet in 50 μ L of elution buffer.
- 13. Vigorously pipet-mix until the beads are uniformly dispersed. Small clumps do not affect performance.
- 14. Incubate at room temperature for **2 minutes** and briefly centrifuge.
- 15. Place the tube on the magnet until the solution is clear, usually ~30 seconds.
- 16. Pipet the eluate (~50 μL) into a new 1.5-mL LoBind[®] tube (purified TCR/BCR PCR1 products).



Purified PCR product can be stored at -20 °C for up to 6 months.

3.3 TCR/BCR PCR2

Summary:

- Prepare TCR/BCR PCR2 mix
- Amplify using TCR/BCR PCR2 program

| Item | | BD Part Number | Preparation and Handling | Storage | | |
|------------------|--|---------------------------------|---|---------|--|--|
| Equili | Equilibrate to room temperature: | | | | | |
| \bigcirc | PCR master mix | 91-1083 | | | | |
| | TCR/BCR universal 91-1205 oligo N2 | | Equilibrate to room temperature 30 minutes before | | | |
| | TCR N2 primer - Mouse OR BCR N2 primer - Mouse | 91-1213 OR 91-1215 | setting up PCR2. Centrifuge briefly. Keep on ice until ready. | −20 °C | | |
| \bigcirc | Nuclease-free water | 650000076 | | | | |
| Obtai | n: | | | | | |
| Purifie | ed TCR/BCR PCR1 product | | | | | |
| Ice bu | Ice bucket | | | | | |
| 0.2-mL PCR tubes | | | | | | |
| Set up | Set up: | | | | | |
| Therm | nocycler with TCR/BCR PCF | R2 program | | | | |

This section describes how to amplify TCR/BCR products through PCR.

1. In the pre-amplification workspace, pipet reagents into a new 1.5-mL LoBind $^{\circledR}$ tube on ice.

TCR and BCR PCR2 mixes

| Сар | Component | For 1 librαry (μL) | For 1 library with 20% overage (µL) | For 4 libraries with 20% overage (µL) | For 8 libraries with 20% overage (µL) |
|-----|---|-----------------------|--|---------------------------------------|---------------------------------------|
| | PCR master mix | 25.0 | 30.0 | 120.0 | 240.0 |
| | TCR/BCR universal oligo N2 | 2.0 | 2.4 | 9.6 | 19.2 |
| | Mouse TCR N2 primer OR Mouse BCR N2 primer ^a | 6.0 | 7.2 | 28.8 | 57.6 |
| | Nuclease-free water | 12.0 | 14.4 | 57.6 | 115.2 |
| | Total | 45.0 | 54.0 | 216.0 | 432.0 |

a. PCR2 mixes for TCR and BCR are made separately.

- 2. Gently vortex mix, briefly centrifuge, and place back on ice.
- 3. Pipet $45 \mu L$ of PCR2 Mix into one 0.2-mL PCR tube for each library.



Bring the TCR PCR2 mix and the BCR PCR2 mix to the post-amplification workspace.

- 4. Pipet 5.0 μL of PCR1 products into 45 μL of PCR2 mix for each library to create the TCR PCR2 Reaction Mix and BCR PCR2 reaction mix, respectively. Total volume of reaction will be $50~\mu L$ for PCR2.
- 5. Gently vortex and briefly centrifuge.
- 6. Run the following PCR program on the thermal cycler.

TCR/BCR PCR2 program

| Step | Cycles | Temperature | Time | |
|-----------|--------|-------------|------------|--|
| Phase I: | 1 | 95 ℃ | 3 minutes | |
| | 15 | 95 ℃ | 30 seconds | |
| | | 70–56 °C | 1 minute | Press Option > Auto |
| | | 72 °C | 1 minute | Delta Starting cycle > "2" Delta > "1 degree" > Done The temperature |
| Phase II: | 8 | 95 ℃ | 30 seconds | |
| | | 55 °C | 1 minute | decreases by 1 °C each cycle, from 70 °C |
| | | 72 °C | 1 minute | to 56 °C. |
| | 1 | 72 °C | 5 minutes | |
| | 1 | 4 °C | ∞ | |



The PCR can run overnight.

3.4 TCR/BCR PCR2 cleanup

Summary:

- TCR/BCR PCR2 cleanup
- Quality check using Qubit Fluorometer

| Item | | BD Part Number | Preparation and Handling | Storage | | | | |
|--|-----------------------------------|----------------|---|---------------------------------|--|--|--|--|
| Equili | Equilibrate to room temperature: | | | | | | | |
| | Elution buffer | 91-1084 | Contain no bairth. | 20.00 | | | | |
| \bigcirc | Nuclease-free water | 65000076 | Centrifuge briefly. | −20 °C | | | | |
| AMPu | re [®] XP magnetic beads | | M 6 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 | | | | | |
| Qubit dsDNA HS Assay Kit Manufacturer's recommendations | | | | | | | | |
| Obtai | n: | | | | | | | |
| TCR/E | 3CR PCR2 product | | | 4 °C | | | | |
| 1.5-m | L DNA LoBind [®] tubes | | | | | | | |
| 0.2-m | L PCR tubes | | | | | | | |
| 0.2-m | 0.2-mL PCR tube magnetic rack | | | | | | | |
| Set u | Set up: | | | | | | | |
| Prepa | re fresh 80% ethyl alcohol | | | Prepare fresh 80% ethyl alcohol | | | | |

This section describes how to perform a single-sided AMPure[®] XP beads cleanup to remove primer dimers from the TCR and BCR PCR2 products. The final product is purified double-stranded DNA.



Perform PCR2 purification in the post-amplification workspace.

1. Make fresh 80% ethyl alcohol for use within 24 hours.



Adjust the volume depending on the number of samples. One sample requires 0.5 mL of 80% ethyl alcohol.

- 2. Bring AMPure[®] XP beads to room temperature and vortex at high speed for **1 minute** until beads are fully resuspended.
- 3. To **50** μ L PCR2 products, pipet **35** μ L of AMPure[®] XP beads (**0.7**x).
- 4. Pipet-mix 10 times and incubate at room temperature for 5 minutes.
- 5. Place the tube on the magnet for 3 minutes. Discard the supernatant.
- Keeping the tube on the magnet, gently add 200 μL of fresh 80% ethyl alcohol into the tube and incubate for 30 seconds.
- 7. Discard the supernatant.
- 8. Repeat steps 6–7 once for a total of two ethyl alcohol washes.
- 9. Keeping the tube on the magnet, use a P20 pipette to remove and discard any residual supernatant from the tube.
- 10. Air-dry the beads at room temperature for 3 minutes.
- 11. Remove the tube from the magnet and resuspend the bead pellet in $50 \mu L$ of elution buffer. Pipet-mix until the beads are fully resuspended.
- 12. Incubate at room temperature for **2 minutes** and briefly centrifuge.
- 13. Place the tube on the magnet until the solution is clear, usually ~30 seconds.
- 14. Pipet the eluate (~50 μ L) into a new 1.5-mL LoBind[®] tube.



Purified PCR product can be stored at -20 °C for up to 6 months.

15. Estimate the concentration by quantifying 2 µL of the TCR/BCR PCR2 library with a Qubit™ Fluorometer using the Qubit™ dsDNA HS Assay Kit. Follow the manufacturer's instructions.

3.5 TCR/BCR RPE

Summary:

- Prepare random priming mix and extension enzyme mix
- Anneal random primers
- Extend random primers

| Item | | BD Part Number | Preparation and Handling | Storage | |
|------------------|---|----------------------|---|---------|--|
| Equili | brate to room temperature: | | | | |
| \bigcirc | TCR/BCR extension buffer | 91-1206 | | | |
| \bigcirc | TCR/BCR extension primers | 91-1208 | | | |
| | dNTP | 650000077 | Equilibrate to room temperature 30 minutes before setting up RPE. Centrifuge briefly. | −20 °C | |
| \bigcirc | Nuclease-free water | 650000076 | | | |
| | Elution buffer | 91-1084 | | | |
| Leave | in freezer until reαdy to use | : | | | |
| 0 | TCR/BCR extension enzyme | 91-1207 | Centrifuge briefly before adding to mix. | −20 °C | |
| Obtai | n: | | | | |
| Purifie | ed TCR/BCR PCR2 product | | | | |
| Ice bu | ıcket | | | | |
| 1.5-m | L DNA LoBind [®] tubes | | | | |
| 0.2-mL PCR tubes | | | | | |
| Set up: | | | | | |
| Therm | Thermocycler with TCR/BCR Denaturation and random priming program | | | | |
| Therm | nocycler with TCR/BCR random | n primer extension p | program | | |

1. Dilute an aliquot of the TCR and BCR PCR2 products with water to 1.0 ng/µL.



If PCR2 concentration is <1 $\,\mathrm{ng}/\mu\mathrm{L}$, increase the volume of PCR2 product needed to ensure 5 $\,\mathrm{ng}$ total concentration and decrease the volume of water in the random primer mix accordingly.

2. In pre-amplification workspace, pipet the following reagents into a new 1.5 mL LoBind $^{\circledR}$ tube:

Random primer mix

| Сар | Component | For 1 librαry (μL) | For 1 library with 20% overage (µL) | For 4 libraries with 20% overage (µL) | For 8 libraries with 20% overage (µL) |
|-----|---------------------------|-----------------------|--|---------------------------------------|---------------------------------------|
| | TCR/BCR extension buffer | 5.0 | 6.0 | 24.0 | 48.0 |
| 0 | TCR/BCR extension primers | 2.5 | 3.0 | 12.0 | 24.0 |
| | Nuclease-free water | Up to 34.0 | Up to 40.8 | Up to 163.2 | Up to 326.4 |
| | Total | 41.5 | 49.8 | 199.2 | 398.4 |

- 3. Pipet-mix the random primer mix and keep at room temperature.
- 4. Pipet 41.5 μ L of random primer mix into one 0.2-mL PCR tube for each library.



Bring the TCR RPE mix and the BCR RPE mix to the post-amplification workspace.

- 5. Add **5** μ L of 1.0 ng/ μ L purified TCR or BCR PCR2 products into each 0.2-mL PCR tube containing random primer mix.
- 6. Total volume of reaction will be $46.5 \mu L$ for random priming.
- 7. Perform denaturation and random priming on thermocycler using the following program:

Program

| Temperature | Time | Cycles |
|-------------|------------|--------|
| 95 ℃ | 5 minutes | |
| 37 ℃ | 5 minutes | 1 |
| 25 ℃ | 15 minutes | |

- 8. Briefly centrifuge the tube and keep at room temperature.
- 9. In pre-amplification workspace, pipet the following reagents into a new 1.5 mL LoBind[®] tube:

Primer extension enzyme mix

| Сар | Component | For 1 librαry (μL) | For 1 library with 20% overage (µL) | For 4 libraries with 20% overage (µL) | For 8 libraries with 20% overage (µL) |
|-----|-----------------------------|-----------------------|--|---------------------------------------|---------------------------------------|
| | dNTP | 2.0 | 2.4 | 9.6 | 19.2 |
| | TCR/BCR Extension Enzyme | 1.5 | 1.8 | 7.2 | 14.4 |
| | Total | 3.5 | 4.2 | 16.8 | 33.6 |

- 10. Gently vortex mix, centrifuge, and place on ice.
- 11. Add $3.5~\mu L$ primer extension enzyme mix to the random priming reaction tube to bring total volume up to **50 \mu L**. Run the following protocol on a thermocycler for extension:

Program

| Temperature | Time | Cycles |
|-------------|------------|--------|
| 25 ℃ | 10 minutes | |
| 37 ℃ | 15 minutes | 1 |
| 45 °C | 10 minutes | |
| 55 ℃ | 10 minutes | |

12. When the PCR program is complete, briefly centrifuge the tubes.

3.6 TCR/BCR RPE cleanup

Summary:

• TCR/BCR RPE cleanup

| Item | | BD Part Number | Preparation and Handling | Storage | | |
|------------|-----------------------------------|----------------|--------------------------------|---------|--|--|
| Equili | Equilibrate to room temperature: | | | | | |
| | Elution buffer | 91-1084 | Contribute by office | −20 °C | | |
| \bigcirc | Nuclease-free water | 65000076 | Centrifuge briefly. | -20 C | | |
| AMPu | re [®] XP magnetic beads | | Manufacturer's recommendations | | | |
| Obtai | n: | | | | | |
| TCR/E | SCR RPE product | | | 4 °C | | |
| 1.5-m | L DNA LoBind [®] tubes | | | | | |
| 0.2-m | L PCR tubes | | | | | |
| 0.2-m | 0.2-mL PCR tube magnetic rack | | | | | |
| Set u | Set up: | | | | | |
| Prepa | Prepare fresh 80% ethyl alcohol | | | | | |



Perform purification in the post-amplification workspace.

1. Make fresh 80% ethyl alcohol for use within 24 hours.



Adjust the volume depending on the number of samples. One sample requires 0.5 mL of 80% ethyl alcohol.

- 2. Bring AMPure® XP beads to room temperature and vortex at high speed for 1 minute until beads are fully resuspended.
- 3. To the **50 \muL** of TCR and BCR RPE products, add **90 \muL** AMPure[®] XP beads (**1.8x**).
- 4. Pipet-mix 10 times and incubate at room temperature for 5 minutes.
- 5. Place the tube on the magnet for 3 minutes.
- 6. Discard the supernatant.
- 7. Keeping the tube on the magnet, gently add $200 \mu L$ of fresh 80% ethyl alcohol into the tube and incubate for 30 seconds.
- 8. Discard the supernatant.
- 9. Repeat steps 7–8 once for a total of two ethyl alcohol washes.
- 10. Keeping the tube on the magnet, use a P20 pipette to remove and discard any residual supernatant from the tube.
- 11. Air-dry the beads at room temperature for 3 minutes.
- 12. Remove tubes from the magnet and add $50 \mu L$ of elution buffer.
- 13. Incubate at room temperature for 2 minutes and briefly centrifuge.
- 14. Place the tube on the magnet until the solution is clear, usually ≤ 30 seconds.
- 15. Pipet the entire eluate (\sim 50 μ L) into a new 1.5-mL LoBind[®] tube separately (purified TCR and BCR RPE products).

3.7 TCR/BCR index PCR

Summary:

- Prepare TCR/BCR index mix
- Amplify using TCR/BCR index program

| Item | | BD Part Number | Preparation and Handling | Storage | | |
|------------|--------------------------------------|----------------|---|---------|--|--|
| Equili | Equilibrate to room temperature: | | | | | |
| | Forward primer 1–8 | Various | | | | |
| | Multiomic reverse primer 1–8 | Various | Equilibrate to room temperature 30 minutes before setting up TCR/BCR Index PCR. | −20 °C | | |
| \bigcirc | Nuclease-free water | 650000076 | Centrifuge briefly. Keep on ice until ready. | | | |
| \bigcirc | PCR master mix | 91-1083 | | | | |
| Obtai | n: | | | | | |
| Purifie | ed TCR/BCR RPE product | | | | | |
| Ice bu | icket | | | | | |
| 1.5-m | 1.5-mL DNA LoBind [®] tubes | | | | | |
| 0.2-m | 0.2-mL PCR tubes | | | | | |
| Set up: | | | | | | |
| Therm | nocycler with TCR/BCR index P | CR program | | | | |

This section describes how to generate TCR/BCR libraries compatible with the Illumina sequencing platform, by adding full-length Illumina sequencing adapters and indices through PCR.

1. In the pre-amplification workspace, pipet reagents into a new 1.5-mL LoBind $^{(8)}$ tube on ice.

TCR/BCR index PCR mix

| Сар | Component | For 1 librαry (μL) | For 1 library with 20% overage (µL) | For 4 libraries with 20% overage (µL) | For 8 libraries with 20% overage (µL) |
|------------|------------------------------|-----------------------|--|---------------------------------------|---------------------------------------|
| \bigcirc | PCR master mix | 25.0 | 30.0 | 120.0 | 240.0 |
| | Forward primer 1–8 | 2.0 | 2.4 | N/A | N/A |
| | Multiomic reverse primer 1–8 | 2.0 | 2.4 | N/A | N/A |
| | Total | 29.0 | 34.8 | 120.0 | 240.0 |

2. Gently vortex mix, briefly centrifuge, and place back on ice.



Bring the TCR/BCR index PCR mix to post-amplification workspace.

- 3. For multiple samples, pipet 25 μ L of index PCR mix into separate 0.2-mL PCR tubes to each sample.
- 4. Add $2 \mu L$ of forward primer and $2 \mu L$ of multiomic reverse primer to each sample.
- 5. Add 21 μ L of TCR/BCR RPE purified products into 29 μ L of TCR/BCR index PCR mix. Total volume of reaction will be $50 \mu L$ for index PCR.
- 6. Gently vortex, and briefly centrifuge.
- 7. Run the following PCR program on the thermal cycler.

TCR/BCR index PCR program

| Step | Cycles | Temperature | Time |
|-----------------|--------|-------------|------------|
| Hot start | 1 | 95 ℃ | 3 minutes |
| Denaturation | | 95 ℃ | 30 seconds |
| Annealing | 10 | 60 °C | 30 seconds |
| Extension | | 72 °C | 30 seconds |
| Final extension | 1 | 72 ℃ | 1 minute |
| Hold | 1 | 4 °C | ∞ |



The PCR can run overnight.

3.8 TCR/BCR index PCR cleanup and quality check

Summary:

- TCR/BCR index PCR cleanup
- Quality check using Qubit Fluorometer and BioAnalyzer/TapeStation

| Item | | BD Part Number | Preparation and Handling | Storage | | | | |
|---|-----------------------------------|----------------|--------------------------------|---------|--|--|--|--|
| Equili | Equilibrate to room temperature: | | | | | | | |
| | Elution buffer | 91-1084 | Contribute briefly | –20 °C | | | | |
| \bigcirc | Nuclease-free water | 650000076 | Centrifuge briefly. | | | | | |
| AMPu | re [®] XP magnetic beads | | | | | | | |
| Qubit | dsDNA HS Assay Kit | | Manufacturer's recommendations | | | | | |
| Agilent BioAnalyzer High Sensitivity Kit OR Agilent TapeStation ScreenTape and Reagents | | | Manufacturer's recommendations | | | | | |
| Obtai | n: | | | | | | | |
| TCR/E | SCR index PCR product | | | 4 °C | | | | |
| 0.2-m | 0.2-mL PCR tubes | | | | | | | |
| 0.2-mL PCR tube magnetic rack | | | | | | | | |
| Set u | Set up: | | | | | | | |
| Prepa | re fresh 80% ethyl alcohol | · | | | | | | |

This section describes how to perform a single-sided AMPure® XP beads cleanup to remove primer dimers from the TCR/BCR index PCR products. The final product is purified double-stranded DNA with full-length Illumina adapter sequences.



Perform index PCR purification in the post-amplification workspace.

1. Make fresh 80% ethyl alcohol for use within 24 hours.



Adjust the volume depending on the number of samples. One sample requires 0.5 mL of 80% ethyl alcohol.

- 2. Bring AMPure® XP beads to room temperature and vortex at high speed for 1 minute until the beads are fully resuspended.
- 3. Briefly centrifuge all the index PCR products.
- 4. Transfer 40 μ L of the TCR and/or BCR index PCR products to a new 0.2-mL PCR tubes.
- 5. Pipet **26 μL** of AMPure[®] XP beads (**0.65x**).
- 6. Pipet-mix 10 times and incubate at room temperature for 5 minutes.
- 7. Place the tube on the strip tube magnet for 3 minutes.
- 8. Discard the supernatant.
- 9. Keeping the tube on the magnet, gently add $200 \mu L$ of fresh 80% ethyl alcohol into the tube and incubate for 30 seconds.
- 10. Discard the supernatant.
- 11. Repeat steps 9–10 for a total of two ethyl alcohol washes.
- 12. Keeping the tube on the magnet, use a P20 pipette to remove and discard the residual supernatant from the tube.
- 13. Air-dry the beads at room temperature for 1 minute.
- 14. Remove the tube from the magnet and resuspend the bead pellet in 50 μ L of elution buffer. Pipet-mix until the beads are fully resuspended.
- 15. Incubate at room temperature for **2 minutes**, and briefly centrifuge.
- 16. Place the tube on the magnet until the solution is clear, usually ~30 seconds.
- 17. Pipet the entire eluate ($\sim 50 \mu L$) into a new 1.5-mL LoBind[®] tube (final sequencing libraries).
- 18. Perform quality control before freezing samples.
 - a. Estimate the concentration by quantifying 2 µL of the final sequencing library with a Qubit™ Fluorometer using the Qubit™ dsDNA HS Kit to obtain an approximate concentration of PCR products to dilute for quantification on an Agilent 2100 Bioanalyzer. Follow the manufacturer's instructions. The expected concentration of the libraries is >1.5 ng/ μ L.

b. Measure the average fragment size of the TCR/BCR library within the size range of 200–1,000 bp by using the Agilent Bioanalyzer with the High Sensitivity Kit for 50–7,000 bp, 5–1,000 pg/ μ L. Follow the manufacturer's instructions.

Figure 2 Sample TapeStation high-sensitivity D5000 trace - Mouse TCR index PCR product





Figure 3 Sample TapeStation high-sensitivity D5000 trace - Mouse BCR index PCR product



Purified PCR product can be stored at –20 °C for up to 6 months.

Sequencing

The sequencing depth for each library is dependent on application. For cell-type clustering, shallow sequencing is sufficient. However, for in-depth analysis, such as comparison across multiple samples, deep sequencing is advised. We recommend meeting the requirement for recursive substitution error correction (RSEC) sequencing depth of ≥6 to reach the threshold of sequencing saturation where most molecules of the library have been recovered, approximately 80%. The RSEC sequencing depth and sequencing saturation are both reported by the analysis pipeline. The actual sequencing reads/cell required to achieve this depth can vary as it is dependent on the gene expression levels, number of cells, and sequencing run quality. The following reads/cell are recommended for WTA mRNA, TCR, and BCR libraries.

Read requirements for libraries

| Library | Read requirement for data analysis | | |
|----------|------------------------------------|--|--|
| WTA mRNA | 10,000–100,000 reads/cell | | |
| TCR | ~5,000 reads/T cell | | |
| BCR | ~5,000 reads/B cell | | |

Required parameters

| Parameter | Requirement |
|------------------|---|
| Platform | Illumina and Element* |
| Paired-end reads | Minimum: Read 1: 65 cycles; Read 2: 150 cycles Recommended: Read 1: 150 cycles; Read 2: 150 cycles Index 1 (i7): 8 cycles; Index 2 (i5): 8 cycles |
| PhiX | Required (3%) |
| Analysis | Refer to the BD [®] Single-Cell Multiomics Bioinformatics Handbook |

^{*} To review Index sequences, see the Appendix (page 73).

Pooling libraries for sequencing

The efficiency of sequencing on Illumina instruments is influenced by many conditions, library size being one of them. The TCR and BCR libraries are ~200–300 bp larger than the WTA mRNA library which will cause them to produce less sequencing data if pooled in a 1:1 ratio with the other libraries. To overcome the difference in sequencing efficiency, more DNA of the TCR and BCR libraries needs to be included in the pool than would be expected when calculating ratios based on read depth. The following tables show examples of different pooling strategies and the expected sequencing outcome, with correction for the size of the TCR and BCR libraries. Validation data indicates a 5X volume correction factor is needed for sequencing TCR and BCR libraries with WTA mRNA libraries.

Example of pooling with correction

In this example, a total of 5,000 enriched T cells were processed. These calculations are using a correction factor of 5 for the TCR library to overcome the differences in sequencing efficiency. The amount of data

needing to be generated (Column D) is based on the cell number (Column B) and expected number of reads per cell (Column C). Based on this example, 150 million reads are needed to achieve the appropriate read depths. Changing the pooling ratios by correcting for the lower TCR sequencing efficiency will help ensure the correct amount of data is generated for each library. This modified pooling scheme, however, does not change the total amount of data needing to be generated, which is 150 million reads.

Pooling for WTA mRNA and TCR libraries

| Α | В | С | D | E | F | G | Н | I | J |
|-----------------|--------------------|---------------------|--------------------------|--|----------------|--------------------------------|-------------------------------------|--------------------|---------------------------------|
| Library type | Number of cells | Expected reads/cell | Reads needed | Pooling ratio before correction | Correction | Reads needed for pooling | Pooling ratio with correction | Sequencing results | Sequencing results (reads/cell) |
| WTA mRNA | 5,000 | 25,000 | 125,000,000 | 83% | N/A | 125,000,000 | 50% | 125,000,000 | 25,000 |
| TCR | 5,000 | 5,000 | 25,000,000 | 17% | 5 ^a | 125,000,000 | 50% | 25,000,000 | 5,000 |
| Total | • | | 150,000,000 ^b | 100% | - | 250,000,000 ^c | 100% | 150,000,000 | - |

a. The 5X correction factor is a recommended starting point and some fine tuning might be required to achieve the optimal library balance.

After sequencing, the total amount of data generated (Column I) as well as the reads/cell for each library (Column J) are as expected (Columns D and C, respectively). The correction for library pooling did not change the amount of data generated (150 million reads) but helped ensure the data was spread out appropriately to each library.

Example of pooling with a mixed population

The following table shows the pooling logic for a mixed population of cells.

Pooling for WTA mRNA, TCR, and BCR libraries

| Α | В | С | D | E | F | G | Н | I | J |
|-----------------|--------------------|---------------------|--------------------------|--|----------------|--------------------------------|-------------------------------------|--------------------|---------------------------------|
| Library type | Number of cells | Expected reads/cell | Reads needed | Pooling ratio before correction | Correction | Reads needed for pooling | Pooling ratio with correction | Sequencing results | Sequencing results (reads/cell) |
| WTA mRNA | 10,000 | 25,000 | 250,000,000 | 88% | N/A | 250,000,000 | 59% | 250,000,000 | 25,000 |
| TCR | 4,000 | 5,000 | 20,000,000 | 7% | 5 ^a | 100,000,000 | 23% | 20,000,000 | 5,000 |
| BCR | 3,000 | 5,000 | 15,000,000 | 5% | 5 ^a | 75,000,000 | 18% | 15,000,000 | 5,000 |
| Total | | | 285,000,000 ^b | 100% | - | 425,000,000 ^c | 100% | 406,000,000 | - |

a. The 5X correction factor is a recommended starting point and some fine tuning might be required to achieve the optimal library balance.

b. Total amount of data to be requested from the sequencing facility plus 3% PhiX.

c. Read total only for pooling purposes.

b. Total amount of data to be requested from the sequencing facility plus 3% PhiX.

c. Read total only for pooling purposes.

Additional considerations

- The 5X volume correction factor for TCR/BCR libraries is a recommended starting place for pooling these libraries. This might need to be adjusted to accommodate different types of flow cells (for example patterned vs non-patterned) or different sequencing platforms.
- · It can be easier to achieve the desired sequencing depth when sequencing multiple TCR or BCR libraries alone since all the libraries are the same size. Pooling will not require a correction and will only be dependent on the number of cells and the reads/cell. This scheme, however, would require 10-15% PhiX, rather than the 3% PhiX when sequencing with the WTA mRNA library.
- Sequencing runs with only WTA libraries require Read 1 of 51 cycles and Read 2 of 71 cycles.
- All libraries derived from the same cartridge can be indexed with the same Illumina indices or reverse index primer from the BD Rhapsody™ reagents. The primary analysis pipeline can differentiate the library types (for example, WTA mRNA vs TCR) based on their structure and sequences. Demultiplexing statistics are reported from the pipeline, but should these statistics be desired prior to running the pipeline, then unique indices will be required for each library.

For additional support with pooling and sequencing, reach out to your local Field Application Specialist (FAS) or scomix@bd.com.

Sequencing recommendations

WTA and TCR/BCR libraries (with or without BD® AbSeg/Sample Tag).

For a NextSeq High or Mid Output run and MiniSeq High or Mid Output run, load the flow cell at a concentration between 1.4–1.8 pM with 3% PhiX for a sequencing run.



If using less than 10,000 reads/cell for the WTA library, increase PhiX percentage to 5–10% to account for lower library diversity.

For NovaSeq:

| Sequencing platform | Cycles | Recommended loading concentration |
|---|----------------------------|--|
| NovaSeq 6000 S Prime (Single Lane) | 2×50, 2×100, 2×150, 2×250* | 180–250 pM (XP workflow) |
| NovαSeq 6000 S Prime (Single Flow Cell) | 2×50, 2×100, 2×150, 2×250* | 350–650 pM (standard workflow) |
| NovaSeq 6000 S1 (Single Lane) | 2×50, 2×100, 2×150* | 180–250 pM (XP workflow) |
| NovαSeq 6000 S1 (Single Flow Cell) | 2×50, 2×100, 2×150* | 350–650 pM (standard workflow) |
| NovαSeq 6000 S2 (Single Flow Cell) | 2×50, 2×100, 2×150* | 350–650 pM (standard workflow) |
| NovaSeq 6000 S4 (Single Lane) | 2×100, 2×150 | 180–250 pM (XP workflow) |
| NovαSeq 6000 S4 (Single Flow Cell) | 2×100, 2×150 | 350–650 pM (standard workflow) |
| NovaSeq X 10B | 2×100, 2×150 | Contact local Field Application Specialist (FAS) |

NovaSeq 100 cycle kit (v1.0 or v1.5) can be used. The 100 cycle kit contains enough reagents for up to 130 cycles.

For other sequencing platforms (e.g. Element AVITI System), follow the manufacturer's sequencing recommendations. Loading concentration may need to be titrated to optimize yield.

Sequencing analysis pipeline

Contact customer support at scomix@bd.com for access to the latest whole transcriptome sequencing analysis pipeline.

Troubleshooting

Library preparation

| Observation | Possible causes | Recommended solutions |
|---|--|--|
| Low yield of indexing PCR. | Input DNA not high enough or cycle number too low. | Repeat indexing PCR with higher cycle number. Alternatively, if RPE-PCR product was diluted before adding to indexing PCR, repeat indexing PCR with less or no dilution. |
| Lower number of reads/cell than expected from mRNA. | 264 bp or ~160 bp products taking reads from mRNA library. | If noise peak is seen in the 264 bp or ~160 bp range, perform a second round of AMPure purification according to 1.8 Additional WTA index PCR cleanup. |

TCR/BCR metrics

| Observation | Possible causes | Recommended solutions |
|--|---|---|
| Low yield of TCR/BCR PCR2, unexpected profile of TCR/BCR | Incorrect components were used. | Ensure BD Rhapsody™ Enhanced Cartridge Reagent Kit V3 (PN 667052) and TSO Next (PN 91-1295) are used for all TCR/BCR Next assays. |
| Index products, or low TCR/BCR pairing efficiency. | Incorrect handling during earlier protocol steps (cDNA synthesis, template switching, denaturation, self-hybridization, and TCR/BCR extension). | Carefully follow all protocol steps in Section 1, especially noting warnings for thermomixer settings, reagent storage temperatures, and incubation timing. |
| | Low viability cells or other challenging samples. | Optimization might be required. Contact your local Field Application Specialist (FAS) or scomix@bd.com. |

Sequencing

| Observation | Possible causes | Recommended solutions |
|--|--|---|
| Over-clustering on the flow cell due to under-estimation of the library. | Inaccurate measurement of the library concentration. | Quantify library according to instructions in protocol. |
| Low sequencing quality. | Insufficient PhiX. | Use the recommended concentration of PhiX with the library to be sequenced. |
| | Suboptimal cluster density, or library denaturation, or both. | See troubleshooting in sequencing platform documentation. |
| One or more libraries are not correctly demultiplexed. | Failure to correctly detect one or more index sequences during sequencing. | Try demultiplexing with a single index. For Illumina sequencers, adjust the mismatch threshold from the default of 1 to allow 2 mismatches. |

| Observation | Possible causes | Recommended solutions |
|--|--|--|
| Failed cluster generation or other sequencing challenges with low-plex pooling. | Sequencing instrument-specific color balance guidelines were not met. For example, the NovaSeq X platform requires signal in the green channel for every cycle. Low % reads in the green channel in a given cycle might not be sufficient. | Consult indexing and pooling guidelines for your sequencing platform. Be sure to consider the final pooling ratio—some libraries might make up a small fraction of the final pool. |

Appendix

Oligonucleotides in BD OMICS-One™ WTA Next Amplification Kit

The following table lists the sequences of all oligonucleotides included in the BD OMICS-One $^{\text{\tiny{M}}}$ WTA Next Amplification Kit (Catalog No. 572620).

| Oligonucleotide | Use | Part/Catalog No. | Sequence (5' – 3') |
|---|--|---------------------|--|
| BD OMICS-One™ Universal Oligo | Forward primer for WTA RPE PCR, Sample Tag PCR1 and PCR2, and BD® AbSeq PCR1 | 51-9025553 | ACA CGA CGC TCT TCC GAT CT |
| BD OMICS-One™ AbSeq Primer | Reverse primer for BD® AbSeq PCR1 | 51-9025468 | CAG ACG TGT GCT CTT CCG ATC T |
| BD OMICS-One™ WTA Extension Primer | Random primer for WTA RPE | 51-9025467 | GGC TCG GAG ATG TGT ATA AGA GAC AG NNNNNNNN |
| BD OMICS-One™ WTA Amplification Primer | Reverse primer for WTA RPE PCR | 51-9025469 | GGC TCG GAG ATG TGT ATA AGA GAC AG |

| Oligonucleotide | Use | Part/Catalog No. | Sequence (5' – 3') |
|---|---|---------------------|---|
| BD OMICS-One™ Library Forward Primer 1 | Forward primer for WTA, Sample Tag, and BD [®] AbSeq | 51-9025472 | AAT GAT ACG GCG ACC ACC GAG ATC TAC AC TATAGCCT ACACTCTTTCCCTACACGACGCTCTTCCGAT*C*T |
| BD OMICS-One™ Library Forward Primer 2 | Index PCR | 51-9025473 | AAT GAT ACG GCG ACC ACC GAG ATC TAC AC ATAGAGGC ACACTCTTTCCCTACACGACGCTCTTCCGAT*C*T |
| BD OMICS-One™ Library Forward Primer 3 | | 51-9025474 | AAT GAT ACG GCG ACC ACC GAG ATC TAC AC CCTATCCT ACACTCTTTCCCTACACGACGCTCTTCCGAT*C*T |
| BD OMICS-One™ Library Forward Primer 4 | | 51-9025475 | AAT GAT ACG GCG ACC ACC GAG ATC TAC AC GGCTCTGA ACACTCTTTCCCTACACGACGCTCTTCCGAT*C*T |
| BD OMICS-One™ Library Forward Primer 5 | | 51-9025476 | AAT GAT ACG GCG ACC ACC GAG ATC TAC AC AGGCGAAG ACACTCTTTCCCTACACGACGCTCTTCCGAT*C*T |
| BD OMICS-One™ Library Forward Primer 6 | | 51-9025477 | AAT GAT ACG GCG ACC ACC GAG ATC TAC AC TAATCTTA ACACTCTTTCCCTACACGACGCTCTTCCGAT*C*T |
| BD OMICS-One™ Library Forward Primer 7 | | 51-9025478 | AAT GAT ACG GCG ACC ACC GAG ATC TAC AC CAGGACGT ACACTCTTTCCCTACACGACGCTCTTCCGAT*C*T |
| BD OMICS-One™ Library Forward Primer 8 | | 51-9025479 | AAT GAT ACG GCG ACC ACC GAG ATC TAC AC GTACTGAC ACACTCTTTCCCTACACGACGCTCTTCCGAT*C*T |

| Oligonucleotide | Use | Part/Catalog No. | Sequence (5' – 3') |
|---|-------------------------------------|---------------------|--|
| BD OMICS-One™ WTA Library Reverse Primer 1 | Reverse primer for WTA Index PCR | 51-9025480 | CAA GCA GAA GAC GGC ATA CGA GAT TACTACGC GTCTCGTGGGCTCGGAGATGTGTATAAGA*G |
| BD OMICS-One™ WTA Library Reverse Primer 2 | | 51-9025600 | CAA GCA GAA GAC GGC ATA CGA GAT AGGCTCCG GTCTCGTGGGCTCGGAGATGTGTATAAGA*G |
| BD OMICS-One™ WTA Library Reverse Primer 3 | | 51-9025482 | CAA GCA GAA GAC GGC ATA CGA GAT GCAGCGTA GTCTCGTGGGCTCGGAGATGTGTATAAGA*G |
| BD OMICS-One™ WTA Library Reverse Primer 4 | | 51-9025483 | CAA GCA GAA GAC GGC ATA CGA GAT CTGCGCAT GTCTCGTGGGCTCGGAGATGTGTATAAGA*G |
| BD OMICS-One™ WTA Library Reverse Primer 5 | | 51-9025484 | CAA GCA GAA GAC GGC ATA CGA GAT GAGCGCTA GTCTCGTGGGCTCGGAGATGTGTATAAGA*G |
| BD OMICS-One™ WTA Library Reverse Primer 6 | | 51-9025485 | CAA GCA GAA GAC GGC ATA CGA GAT CGCTCAGT GTCTCGTGGGCTCGGAGATGTGTATAAGA*G |
| BD OMICS-One™ WTA Library Reverse Primer 7 | | 51-9025486 | CAA GCA GAA GAC GGC ATA CGA GAT GTCTTAGG GTCTCGTGGGCTCGGAGATGTGTATAAGA*G |
| BD OMICS-One™ WTA Library Reverse Primer 8 | | 51-9025487 | CAA GCA GAA GAC GGC ATA CGA GAT ACTGATCG GTCTCGTGGGCTCGGAGATGTGTATAAGA*G |

| Oligonucleotide | Use | Part/Catalog No. | Sequence (5' – 3') |
|---|--|---------------------|--|
| BD OMICS-One™ Multiomic Library Reverse Primer 1 | Reverse primer for Sample Tag and BD [®] AbSeq Index PCR | 51-9025489 | CAA GCA GAA GAC GGC ATA CGA GAT TACTACGC GTGACTGGAGTTCAGACGTGTGCTCTTCCGATC*T |
| BD OMICS-One™ Multiomic Library Reverse Primer 2 | | 51-9025490 | CAA GCA GAA GAC GGC ATA CGA GAT AGGCTCCG GTGACTGGAGTTCAGACGTGTGCTCTTCCGATC*T |
| BD OMICS-One™ Multiomic Library Reverse Primer 3 | | 51-9025492 | CAA GCA GAA GAC GGC ATA CGA GAT GCAGCGTA GTGACTGGAGTTCAGACGTGTGCTCTTCCGATC*T |
| BD OMICS-One™ Multiomic Library Reverse Primer 4 | | 51-9025493 | CAA GCA GAA GAC GGC ATA CGA GAT CTGCGCAT GTGACTGGAGTTCAGACGTGTGCTCTTCCGATC*T |
| BD OMICS-One™ Multiomic Library Reverse Primer 5 | | 51-9025494 | CAA GCA GAA GAC GGC ATA CGA GAT GAGCGCTA GTGACTGGAGTTCAGACGTGTGCTCTTCCGATC*T |
| BD OMICS-One™ Multiomic Library Reverse Primer 6 | | 51-9025496 | CAA GCA GAA GAC GGC ATA CGA GAT CGCTCAGT GTGACTGGAGTTCAGACGTGTGCTCTTCCGATC*T |
| BD OMICS-One™ Multiomic Library Reverse Primer 7 | | 51-9025497 | CAA GCA GAA GAC GGC ATA CGA GAT GTCTTAGG GTGACTGGAGTTCAGACGTGTGCTCTTCCGATC*T |
| BD OMICS-One™ Multiomic Library Reverse Primer 8 | | 51-9025498 | CAA GCA GAA GAC GGC ATA CGA GAT ACTGATCG GTGACTGGAGTTCAGACGTGTGCTCTTCCGATC*T |

| Forward Index Name | i5 bases for sample sheet | i5 bases for sample sheet |
|--|------------------------------------|--|
| | NovaSeq, MiSeq, HiSeq 2000/2500 | iSeq, MiniSeq, NextSeq, HiSeq 3000/4000 |
| BD OMICS-One™ Library Forward Primer 1 | TATAGCCT | AGGCTATA |
| BD OMICS-One™ Library Forward Primer 2 | ATAGAGGC | GCCTCTAT |
| BD OMICS-One™ Library Forward Primer 3 | ССТАТССТ | AGGATAGG |
| BD OMICS-One™ Library Forward Primer 4 | GGCTCTGA | TCAGAGCC |
| BD OMICS-One™ Library Forward Primer 5 | AGGCGAAG | СТТССССТ |
| BD OMICS-One™ Library Forward Primer 6 | TAATCTTA | TAAGATTA |
| BD OMICS-One™ Library Forward Primer 7 | CAGGACGT | ACGTCCTG |
| BD OMICS-One™ Library Forward Primer 8 | GTACTGAC | GTCAGTAC |

| Reverse Index Name | Bases in adapter | i7 bases for sample sheet |
|--|------------------|---------------------------|
| BD OMICS-One™ WTA Library Reverse Primer 1 BD OMICS-One™ Multiomic Library Reverse Primer 1 | TACTACGC | GCGTAGTA |
| BD OMICS-One™ WTA Library Reverse Primer 2 BD OMICS-One™ Multiomic Library Reverse Primer 2 | AGGCTCCG | CGGAGCCT |
| BD OMICS-One™ WTA Library Reverse Primer 3 BD OMICS-One™ Multiomic Library Reverse Primer 3 | GCAGCGTA | TACGCTGC |
| BD OMICS-One™ WTA Library Reverse Primer 4 BD OMICS-One™ Multiomic Library Reverse Primer 4 | CTGCGCAT | ATGCGCAG |
| BD OMICS-One™ WTA Library Reverse Primer 5 BD OMICS-One™ Multiomic Library Reverse Primer 5 | GAGCGCTA | TAGCGCTC |
| BD OMICS-One™ WTA Library Reverse Primer 6 BD OMICS-One™ Multiomic Library Reverse Primer 6 | CGCTCAGT | ACTGAGCG |
| BD OMICS-One™ WTA Library Reverse Primer 7 BD OMICS-One™ Multiomic Library Reverse Primer 7 | GTCTTAGG | CCTAAGAC |
| BD OMICS-One™ WTA Library Reverse Primer 8 BD OMICS-One™ Multiomic Library Reverse Primer 8 | ACTGATCG | CGATCAGT |

Mouse T cell PCR1 primers

| Primer name | Primer sequence (5'–3') |
|-----------------------|-----------------------------|
| Ms_TRAC_N1 | TTTTCGGCACATTGATTTGGGAG |
| Ms_TRBC_N1 | CTCAGGCAGTAGCTATAATTGCT |
| Ms_TRDC_N1 | CAATCTTCTTGGATGATCTGAGACT |
| Ms_TRGC1- TRGC2_N1 | GGAAAGAACTTTTCAAGGAGACAAAGG |

Mouse T cell PCR2 primers

| Primer name | Primer sequence (5'-3') |
|-----------------------|----------------------------|
| Ms_TRAC_N2 | AGGTTCTGGGTTCTGGATGT |
| Ms_TRBC_N2 | CAATCTCTGCTTTTGATGGCTC |
| Ms_TRDC_N2 | GTAGAAATCTTTCACCAGACAAGC |
| Ms_TRGC1- TRGC2_N2 | TTGGGGGAAATGTCTGCA |
| Ms_TRGC4_N2 | ATAGTAGGCTTGGGAGAAAAGTCTGA |

Mouse B cell PCR1 primers

| Primer name | Primer sequence (5'–3') |
|-------------------------|--------------------------|
| Ms_IGHA_N1 | AACTGGCTGCTCATGGTGTA |
| Ms_IGHD_N1 | AAGTGTGGTTGAGGTTCAGTTCTG |
| Ms_IGHE_N1 | GAAGTTCACAGTGCTCATGTTC |
| Ms_IGHG1_N1 | CAGAGTGTAGAGGTCAGACT |
| Ms_IGHG2A- IGHG2C_N1 | TCGAGGTTACAGTCACTGAG |
| Ms_IGHG2B_N1 | GATCCAGAGTTCCAAGTCACAG |
| Ms_IGHG3_N1 | TACGTTGCAGATGACAGTCT |
| Ms_IGHM_N1 | TGGATGACTTCAGTGTTGTTCTG |
| Ms_IGKC_N1 | TGTAGGTGCTGTCTTTGCTG |
| Ms_IGLC1_N1 | CTGTAACTGCTATGCCTTTCCC |
| Ms_IGLC2-IGLC3_ N1 | TTGGTGGGATTTGAAGTGTCC |

Mouse B cell PCR2 primers

| Primer name | Primer sequence (5'–3') |
|-------------------------|--------------------------|
| Ms_IGHA_N2 | TGTCAGTGGGTAGATGGTGG |
| Ms_IGHD_N2 | CTGACTTCCAATTACTAAACAGCC |
| Ms_IGHE_N2 | TAGAGCTGAGGGTTCCTGATAG |
| Ms_IGHG1_N2 | CAGTGGATAGACAGATGGGGGT |
| Ms_IGHG2A- IGHG2C_N2 | ATGGGGCTGTTGTTTTGG |
| Ms_IGHG2B_N2 | GTGGATAGACTGATGGGGGTGTT |
| Ms_IGHG3_N2 | AGGGAAGTAGCCTTTGACAAG |
| Ms_IGHM_N2 | GACATTTGGGAAGGACTGACTC |
| Ms_IGKC_N2 | AGATGTTAACTGCTCACTGGATG |
| Ms_IGLC1_N2 | GTTAGTCTCGAGCTCTTCAGA |
| Ms_IGLC2-IGLC3_N2 | CAGTGTGGCTTTGTTTTCCT |

Contact Information

Becton, Dickinson and Company BD Biosciences 155 North McCarthy Boulevard Milpitas, California 95035 USA

bdbiosciences.com scomix@bd.com