

# Single-Cell Capture and cDNA Synthesis with the BD Rhapsody™ Express Single-Cell Analysis System

For Research Use Only

Doc ID: 210967 Rev. 3.0  
23-22952(01)  
2022-01



---

Becton, Dickinson and Company  
BD Biosciences  
2350 Qume Drive  
San Jose, California 95131 USA

bdbiosciences.com  
scomix@bdscomix.bd.com

## Copyrights

No part of this publication may be reproduced, transmitted, transcribed, stored in retrieval systems, or translated into any language or computer language, in any form or by any means: electronic, mechanical, magnetic, optical, chemical, manual, or otherwise, without prior written permission from BD.

The information in this document is subject to change without notice. BD reserves the right to change its products and services at any time to incorporate the latest technological developments. Although this guide has been prepared with every precaution to ensure accuracy, BD assumes no liability for any errors or omissions, nor for any damages resulting from the application or use of this information. BD welcomes customer input on corrections and suggestions for improvement.

## Trademarks

BD, the BD Logo and Rhapsody are trademarks of Becton, Dickinson and Company or its affiliates. All other trademarks are the property of their respective owners. © 2022 BD. All rights reserved.

## Regulatory Information

For Research Use Only. Not for use in diagnostic or therapeutic procedures.

## History

Revision	Date	Change made
210967 Rev. 1.0	2018-07	Initial release
23-22952-00 (210967 Rev. 2.0)	2020-07	Assigned new 23 document part number, updated catalog number for reverse transcriptase, recommended thermomixer instead of heat block for Exonuclease I inactivation
23-22952(01) (210967 Rev. 3.0)	2022-01	Added Enhanced Cell Capture Beads and part numbers.

# Contents

- Introduction ..... 5
  - Required materials ..... 5
  - Best practices ..... 5
  - Before you begin ..... 6
- Priming and treating the BD Rhapsody™ Cartridge ..... 6
- Counting and preparing a single-cell suspension for cartridge loading ..... 7
- Loading cells in the cartridge ..... 8
- Preparing Cell Capture Beads ..... 9
- Loading and washing Cell Capture Beads ..... 9
- Lysing cells ..... 10
- Retrieving Cell Capture Beads ..... 11
- Washing Cell Capture Beads ..... 11
- Performing reverse transcription ..... 12
- Treating the sample with Exonuclease I ..... 13
- Troubleshooting ..... 14



# Introduction

This protocol describes cell loading in the BD Rhapsody™ Cartridge and single-cell capture with the BD Rhapsody™ Express Single-Cell Analysis system.

For complete instrument procedures, safety information, and workflow, see the *BD Rhapsody™ Express Single-Cell Analysis System Instrument User Guide* (Doc ID: 214063).

## Required materials

For a complete list of materials, see instrument user guide.

- BD Rhapsody™ Cartridge Reagent Kit (Cat. No. 633731) or BD Rhapsody™ Enhanced Cartridge Reagent Kit (Cat. No. 664887)
- BD Rhapsody™ Cartridge Kit (Cat. No. 633733)
- BD Rhapsody™ cDNA Kit (Cat. No. 633773)
- Falcon® Tube with Cell Strainer Cap (Corning Cat. No. 352235)
- BD Rhapsody™ P1200M pipette (Cat. No. 633704)
- BD Rhapsody™ P5000M pipette (Cat. No. 633705)
- Large magnetic separation stand (V&P Scientific, Inc. Cat. No. VP 772FB-1)
- 15-mL tube adapter (V&P Scientific Cat. No. VP 772FB-1A)
- 6-Tube Magnetic Separation Rack for 1.5-mL tubes (New England Biolabs Cat. No. S1506S)

## Best practices

- Always use low-retention filtered pipette tips and LoBind tubes.
- Perform single-cell capture and cDNA synthesis in a pre-amplification workspace.
- Prepare cells as close to cell loading as possible. Keep the other reagents, including Sample Buffer (Cat. No. 650000062), on ice, unless instructed otherwise.
- Change pipette tips before every pipetting step.
- To ensure an air-tight seal with the BD Rhapsody™ P1200M (Cat. No. 633704) and P5000M (Cat. No. 633705) pipettes, hold the pipette with one hand and slightly twist the pipette to firmly seat a pipette tip on the pipette shaft.

## Before you begin

- Thaw reagents (not enzymes) in the BD Rhapsody™ cDNA Kit (Cat. No. 633773) at room temperature (15°C to 25°C), and then place on ice. Keep enzymes at –25°C to –15°C.
- Place on ice the following components of the BD Rhapsody™ Cartridge Reagent Kit:
  - Sample Buffer (Cat. No. 650000062)
  - 1 M DTT (Cat. No. 650000063)
  - Lysis Buffer (Cat. No. 650000064)
  - Cell Capture Beads (Cat. No. 650000089) or BD Rhapsody™ Enhanced Cell Capture Beads (Cat. No. 700027881)
- Ensure that the SmartBlock™ Thermoblock 1.5 mL or equivalent is installed on the thermomixer and is set to 37°C for 20 minutes.
- Set a additional thermomixer to 80°C.
- Prepare a single-cell suspension. See *Preparing Single-Cell Suspensions Protocol* (Doc ID: 210964).
- If your biological sample contains red blood cell contamination, red blood cell lysis is required. See *Preparing Single-Cell Suspensions Protocol* (Doc ID: 210964).

## Priming and treating the BD Rhapsody™ Cartridge

Prime and treat the BD Rhapsody™ Cartridge (Cat. No. 400000847). For detailed instructions, see the instrument user guide.

Express instrument slider	Position
Front	Waste
Side	0

Step no.	Material to load	Volume (µL)	P1200M pipette mode	Incubation at room temperature
1	100% ethyl alcohol	700	Prime/Treat	—
2	Air	700	Prime/Treat	—
3	Room temp. Cartridge Wash Buffer 1 (Cat. No. 650000060)	700	Prime/Treat	1 min
4	Air	700	Prime/Treat	—
5	Room temp. Cartridge Wash Buffer 1 (Cat. No. 650000060)	700	Prime/Treat	10 min

Step no.	Material to load	Volume (µL)	P1200M pipette mode	Incubation at room temperature
6	Air	700	Prime/Treat	—
7	Room temp. Cartridge Wash Buffer 2 (Cat. No. 650000061)	700	Prime/Treat	≤4 hr

## Counting and preparing a single-cell suspension for cartridge loading

For detailed instructions, see the instrument user guide.

- 1 Treat the cells with viability stain(s), and count. Order of accurate counting:
  - Manual counting with fluorescence
  - Automated counting with fluorescence
  - Automated counting with Trypan Blue Stain and brightfield
  - Manual counting with Trypan Blue Stain and brightfield
- 2 Determine the desired number of cells to capture in the BD Rhapsody™ Cartridge. See the instrument user guide for a table containing estimated multiplet rates based on the number of captured cells on retrieved Cell Capture Beads.
- 3 Determine the pooling ratio of samples to load onto the BD Rhapsody™ Cartridge. For example, if two samples were labeled using the BD Rhapsody™ Single-Cell Multiplexing Kit (Cat. No. 633781), and the samples will be pooled in equal proportion, the pooling ratio for each sample is 0.5. If only one sample is used, the pooling ratio is 1.
- 4 Calculate the volume,  $V$ , for each sample needed to prepare the pooled single-cell suspension:

$$V = N \times P \times 1.36 / C$$

where:

$V$  = volume of sample needed (µL)

$N$  = desired number of captured cells in cartridge

$P$  = pooling ratio

$C$  = total cell concentration (cells/µL)

### Example

On a BD Rhapsody™ Cartridge, you want to capture 10,000 cells that are pooled equally of Sample A and Sample B.

$N$  = desired number of captured cells in cartridge = 10,000

$P_A$  = sample A pooling ratio = 0.5

$P_B$  = sample B pooling ratio = 0.5

$C_A$  = sample A total cell concentration = 200 cells/ $\mu$ L  
 $C_B$  = sample B total cell concentration = 400 cells/ $\mu$ L

$$\text{Volume of sample A needed} = 10,000 \text{ cells} \times 0.5 \times 1.36 / 200 \text{ cells}/\mu\text{L} = 34 \mu\text{L}$$

$$\text{Volume of sample B needed} = 10,000 \text{ cells} \times 0.5 \times 1.36 / 400 \text{ cells}/\mu\text{L} = 17 \mu\text{L}$$

- 5 Calculate the sum of all of the sample volumes,  $V_n$ , to be used in the cell suspension. Using the example in [step 4](#):

$$V_n = 34 \mu\text{L} + 17 \mu\text{L} = 51 \mu\text{L}$$

- 6 Calculate the volume of cold Sample Buffer,  $B$ , that is needed to bring the final volume of cell suspension to 650  $\mu$ L. Using the example in [step 5](#):

$$B = 650 \mu\text{L} - 51 \mu\text{L} = 599 \mu\text{L}$$

**NOTE** For low-abundance samples, the final cell suspension can be prepared in 610  $\mu$ L of cold Sample Buffer.

- 7 According to the calculations in steps 3–6, prepare the cell suspension in cold Sample Buffer (Cat. No. 650000062) in a new 1.5-mL LoBind tube.

Ensure the stock solution is well resuspended by gentle pipet-mixing before pooling.

- 8 If the samples were not filtered before counting the cells, filter through a Falcon<sup>®</sup> tube with a Cell Strainer Cap (Corning Cat. No. 352235).

## Loading cells in the cartridge

- 1 Load the cartridge with materials listed using the P1200M pipette:

Material to load	Volume ( $\mu$ L)	Pipette mode
Air	700	Prime/Treat
<ul style="list-style-type: none"> <li>Set P1200M pipette to Cell Load mode.</li> <li>Pipet-mix the cell suspension with a manual P1000 pipette.</li> </ul>		
Cell suspension	575	Cell Load <sup>a</sup>

- a. Press button once to aspirate 40  $\mu$ L air, and then immerse tip in cell suspension. Press button again to aspirate 575  $\mu$ L of cold cell suspension. Dispense 615  $\mu$ L of air and cell suspension.



**NOTE** Air bubbles that might appear at the inlet or outlet of the cartridge do not affect cartridge performance.

- 2 If necessary, wipe condensation from the top cartridge surface for optimal scanning.
- 3 Incubate at room temperature (15°C to 25°C) for 15 minutes.

During the 15-minute incubation, prepare the Cell Capture Beads (Cat. No. 650000089) or BD Rhapsody™ Enhanced Cell Capture Beads (Cat. No. 700027881). See [Preparing Cell Capture Beads](#) in the following section.

## Preparing Cell Capture Beads

Keep the Cell Capture Beads on ice before use.

For maximum recovery, do not vortex samples containing Cell Capture Beads. Gently mix the suspensions with Cell Capture Beads by pipette only. Use low-retention pipette tips and LoBind tubes. Keep the beads cold, and pipet-mix only.

- 1 Place the Cell Capture Bead tube on the magnet for 1 minute, and remove the storage buffer.
- 2 Remove the tube from the magnet, and pipet 750 µL of cold Sample Buffer (Cat. No. 650000062) into the tube.
- 3 Pipet-mix, and place on ice.

## Loading and washing Cell Capture Beads

- 1 Set the P1200M pipette to Prime/Treat mode.
- 2 Load the cartridge with the following materials using the P1200M pipette:

Material to load	Volume (µL)	Pipette mode
Air	700	Prime/Treat
<ul style="list-style-type: none"><li>• Set P1200M pipette to <b>Bead Load</b> mode.</li><li>• Use a manual P1000 to gently pipet-mix the beads in cold Sample Buffer (Cat. No. 650000062). Immediately load.</li></ul>		
Cell Capture Beads	630	Bead Load

- 3 Incubate the cartridge at room temperature (15°C to 25°C) for 3 minutes.
- 4 Place the cartridge on the plate shaker plate adapter.
- 5 Shake the cartridge at room temperature (15°C to 25°C) for 15 seconds.

Following is a list of suggested shakers/mixers with their corresponding settings:

- Eppendorf ThermoMixer® C: 1,000 rpm
- Eppendorf MixMate®: 1,000 rpm
- MicroPlate Genie®: 1,600 rpm. Set external timer to 15 seconds

- 6 Blot the outlet drip with a lint-free wipe.
- 7 Return the cartridge to the Express instrument, and wait 30 seconds.
- 8 Set the P1200M pipette to **Wash** mode.
- 9 Load the cartridge with the following materials using the P1200M pipette:

Material to load	Volume (µL)	Pipette mode <sup>a</sup>
Air	700	Wash
Cold Sample Buffer (Cat. No. 650000062)	700	Wash
Air	700	Wash
Cold Sample Buffer (Cat. No. 650000062)	700	Wash

a. Press button once to aspirate 720 µL air or reagent. Insert the tip into the cartridge, and press button once to dispense 700 µL air or liquid. Remove pipette tip, and press button once to dispense remaining 20 µL of air or liquid.

## Lysing cells

Avoid bubbles.

- 1 Add 75.0 µL of 1 M DTT (Cat. No. 650000063) to one 15-mL Lysis Buffer bottle (Cat. No. 650000064). Check the **Add DTT** box.

Use the Lysis Buffer with DTT ≤24 hours, and then discard.

- 2 Briefly vortex the lysis mix, and place on ice.
- 3 Move the left slider to **LYSIS** on the Express instrument.
- 4 Set the P1200M pipette to **Lysis** mode.
- 5 Load the cartridge with the following material using the P1200M pipette:

Material to load	Volume (µL)	Pipette mode
Lysis Buffer with DTT	550	Lysis

- 6 Incubate at room temperature (15°C to 25°C) for 2 minutes.  
Maintain the recommended lysis time for best performance.

## Retrieving Cell Capture Beads

- 1 Place the 5-mL LoBind tube in the Express instrument drawer.
- 2 Ensure that the P5000M pipette is set to **Retrieval** mode.
- 3 Move the front slider to **BEADS** on the Express instrument.
- 4 Move the left slider to **RETRIEVAL**.
- 5 Leave the Retrieval magnet in the down position for 30 seconds.
- 6 Aspirate 5,000 µL of Lysis Buffer with DTT using the P5000M pipette.
- 7 Press down on the P5000M pipette to seal against the gasket.
- 8 Move the left slider to the middle position (0), and *immediately* load 4,950 µL of Lysis Buffer with DTT.
- 9 Remove the pipette from the gasket, and purge the tip.
- 10 Move the front slider to **OPEN**, and place the 5-mL LoBind tube on the large magnet with the 15-mL tube adapter (V&P Scientific Cat. No. VP 772FB-1A) for 1 minute.
- 11 Immediately proceed to [Washing Cell Capture Beads](#) in the following section.
- 12 Appropriately dispose of the cartridge, Waste Collection Container, and Lysis Buffer with DTT.



**Biological hazard.** All surfaces that come in contact with biological specimens can transmit potentially fatal disease. Use universal precautions when cleaning surfaces. Wear suitable protective clothing, eyewear, and gloves.

- 13 Clean the Express instrument with 10% bleach or 70% ethyl alcohol.

## Washing Cell Capture Beads

- 1 After the 1-minute incubation, leaving the 5-mL tube containing retrieved Cell Capture Beads on the large magnet, remove all but ~1 mL of the supernatant without disturbing the beads.
- 2 Remove the tube from the magnet. Gently pipet-mix the beads and transfer them to a new 1.5-mL LoBind tube.
- 3 If there are still beads left in the 5-mL tube, add 0.5 mL of Lysis Buffer with DTT, rinse the 5-mL tube, and transfer to the 1.5-mL LoBind tube from the previous step.

- 4 Place the tube on the 1.5-mL magnet for  $\leq 2$  minutes. Remove the supernatant.  
Avoid leaving Lysis Buffer or bubbles in the tube. Lysis Buffer might cause the reverse transcription reaction to fail.
- 5 Remove the tube from the magnet, and pipet 1.0 mL of cold Bead Wash Buffer (Cat. No. 650000065) into the tube. Pipet-mix.
- 6 Place the tube on the 1.5-mL magnet for  $\leq 2$  minutes. Remove the supernatant.
- 7 Remove the tube from the magnet, and pipet 1.0 mL of cold Bead Wash Buffer (Cat. No. 650000065) into the tube. Pipet-mix, and place on ice.

Start reverse transcription  $\leq 30$  minutes after washing the retrieved Cell Capture Beads with Bead Wash Buffer.

## Performing reverse transcription

- 1 Ensure that the SmartBlock™ Thermoblock for ThermoMixer® C is at 37°C.
- 2 In the pre-amplification workspace, pipet the following reagents into a new 1.5-mL LoBind tube on ice:

### cDNA mix

Component	For 1 library (μL)	For 1 library + 20% overage (μL)
RT Buffer (Cat. No. 650000067)	40.0	48.0
dNTP (Cat. No. 650000077)	20.0	24.0
RT 0.1 M DTT (Cat. No. 650000068)	10.0	12.0
Bead RT/PCR Enhancer (Cat. No. 91-1082)	12.0	14.4
RNase Inhibitor (Cat. No. 650000078)	10.0	12.0
Reverse Transcriptase (Cat. No. 700026321)	10.0	12.0
Nuclease-Free Water (Cat. No. 650000076)	98.0	117.6
<b>Total</b>	<b>200.0</b>	<b>240.0</b>

- 3 Gently vortex mix, briefly centrifuge, and place back on ice.
- 4 Place the tube of washed Cell Capture Beads on the 1.5-mL tube magnet for  $\leq 2$  minutes. Remove the supernatant.
- 5 Remove the tube from the magnet and pipet 200 μL of cDNA mix into the beads. Pipet-mix.  
Keep the prepared cDNA mix with beads on ice until the suspension is transferred in the next step.

- 6 Transfer the bead suspension to a new 1.5-mL LoBind tube.
- 7 Incubate the bead suspension on the SmartBlock™ Thermoblock for ThermoMixer® C at 1,200 rpm and 37°C for 20 minutes.  
  
Shaking is critical for this incubation.
- 8 Place on ice.

## Treating the sample with Exonuclease I

- 1 Set one thermomixer to 37°C and a second thermomixer to 80°C.

**NOTE** Exonuclease I inactivation temperatures above 80°C can result in the loss of AbSeq molecules, thus a heat block should not be used for this step. If only one thermomixer is available, allow it to equilibrate to 80°C before starting the inactivation incubation.

- 2 In the pre-amplification workspace, pipet the following reagents into a new 1.5-mL LoBind tube on ice:

### Exonuclease I mix

Component	For 1 library (µL)	For 1 library + 20% overage (µL)
10X Exonuclease I Buffer (Cat. No. 650000071)	20.0	24.0
Exonuclease I (Cat. No. 650000072)	10.0	12.0
Nuclease-Free Water (Cat. No. 650000076)	170.0	204.0
<b>Total</b>	<b>200.0</b>	<b>240.0</b>

- 3 Gently vortex mix, briefly centrifuge, and place back on ice.
- 4 Place the tube of Cell Capture Beads with cDNA mix on the 1.5-mL tube magnet for ≤2 minutes. Remove the supernatant.
- 5 Remove the tube from the magnet, and pipet 200 µL of Exonuclease I mix into the tube. Pipet-mix.
- 6 Incubate the bead suspension on the thermomixer at 1,200 rpm and 37°C for 30 minutes.  
  
If only one thermomixer is available, allow it to equilibrate to 80°C before starting the inactivation incubation. Place the samples on ice until that temperature is reached.
- 7 Incubate the bead suspension on the thermomixer (no shaking) at 80°C for 20 minutes.  
  
Do not exceed this inactivation temperature and incubation time.
- 8 Place the tube on ice for ~1 minute.

- 9 Place the tube on the magnet for  $\leq 1$  minute until clear. Remove the supernatant.
- 10 Remove the tube from the magnet, and pipet 200  $\mu\text{L}$  of cold Bead Resuspension Buffer (Cat. No. 650000066) into the tube. Pipet-mix.

**Stopping point:** Exonuclease I-treated beads can be stored at 2°C to 8°C for up to 3 months.

- 11 Proceed to library preparation. See the *Single-Cell Analysis Workflow with BD Rhapsody™ Systems* (Doc ID: 220524).

## Troubleshooting

For additional troubleshooting on scanning or cartridge loading, see the troubleshooting section in the instrument user guide.

For technical support, contact [scmix@bdscomix.bd.com](mailto:scmix@bdscomix.bd.com).

Observation	Possible causes	Recommended solutions
No pellet after centrifuging cells or very few cells	Rare or dilute sample	After each centrifugation step, leave 50 $\mu\text{L}$ of supernatant.