

Intracellular CITE-seq using BD[®] AbSeq Antibody-Oligos

Workflow presentation



Redefine CITE-seq: Profile intracellular proteins, surface proteins and RNA in one single-cell experiment

Single-cell multiomics

- Simultaneously detect single-cell transcriptome, surface and intracellular proteome profiles
- Yields high WTA sensitivity* for RNA analyses
- Faithfully recover expected cell subsets in your sample

High-parameter proteomics profiling

- Profile up to 100 protein markers including surface and intracellular proteins in a single experiment

Sample multiplexing enabled

- Multiplex up to 12 samples using the BD® Single-Cell Multiplexing Kit (SMK)

Validated assay

- Fully validated on the BD Rhapsody™ System
- Each IC BD® AbSeq Ab-Oligo is validated against flow cytometry controls
- Comparable IC protein performance vs flow cytometry

Simple and free bioinformatics analysis

- Analyze data with the intuitive BD Rhapsody™ Sequence Analysis Pipeline
- Same pipeline for analysis of intracellular and surface proteins

Efficient and flexible workflow

- ~2 hours in addition to the regular BD Rhapsody™ System workflow
- Stopping point for up to 24 hours during the workflow without compromising assay sensitivity

Required reagents for the IC CITE-seq Assay using BD[®] AbSeq Ab-Oligos

Required reagents for CITE-seq with BD[®] AbSeq Ab-Oligos

BD Rhapsody™ Enhanced Cartridge Reagent Kit

Cat. No. 664887, for four cartridge lane runs

BD Rhapsody™ Cartridge Kit

Cat. No. 633733, four cartridge lanes

OR

BD Rhapsody™ 8-Lane Cartridge

Cat. No. 666262, eight cartridge lanes

BD Rhapsody™ cDNA Kit

Cat. No. 633773, for four cartridge lane runs

BD Rhapsody™ WTA Amplification Kit

Cat. No. 633801, for four cartridge lane runs

Order the single-vial surface BD[®] AbSeq Ab-Oligos separately to suit experimental needs

Additional reagents needed for IC CITE-seq Assay

BD Rhapsody™ Intracellular AbSeq Buffer Kit

Cat. No. 570742, for four IC staining tests

Base buffer, 20 mL

Perm buffer, 12 mL

0.1 M DTT, 20 µL

Proteinase K, 210 µL

Nuclease-free water, 20 mL

BD[®] RNase Inhibitor

Cat. No. 570751, for four IC staining tests

Single-vial reagent, 900 µL

OR

RNase Inhibitor (New England Biolabs, M0314L)

Using RI from other vendors could impact assay performance

BD[®] AbSeq Enhancer

Cat. No. 570750, for four IC staining tests*

BD[®] AbSeq Enhancer 1, 50 µL

BD[®] AbSeq Enhancer 2, 50 µL

BD[®] AbSeq Enhancer 3, 50 µL

BD[®] OMICS-Guard Sample Preservation Buffer Kit

Cat. No. 570908, for 12 million cells

12 vials/kit, 1 test/vial (1 mL)

OR

BD[®] OMICS-Guard Sample Preservation Buffer

Cat. No. 570911

50-mL bottle

Need to aliquot in sterile conditions

Vybrant™ DyeCycle™ Green Stain (Thermo Fisher Scientific, V35004)

Order the single-vial IC BD[®] AbSeq Ab-Oligos separately to suit experimental needs

Workflow

IC CITE-seq Assay workflow

Surface CITE-seq Assay workflow



IC CITE-seq Assay workflow



Intracellular BD[®] AbSeq Ab-Oligos and supporting products



40 min

40 min

80 min

Hands-on time




120 min

130 min

90 min

Total time

IC CITE-seq Assay chemistry




Surface AbSeq Staining

- 1 **Surface staining:** Surface AbSeq/SMK bind to surface proteins


IC CITE-seq Assay chemistry

SAFE STOPPING POINT:

User will have an option to store the cells for up to 24 h in BD® OMICS-Guard Buffer before moving to permeabilization



Surface AbSeq Staining




Preservation

- 1 **Surface staining:** Surface AbSeq/SMK bind to surface proteins
- 2 **Preservation:** Preserves cell status


IC CITE-seq Assay chemistry

SAFE STOPPING POINT:


User will have an option to store the cells for up to 24 h in BD® OMICS-Guard Buffer before moving to permeabilization



Surface AbSeq Staining



Preservation




Permeabilization

- 1 **Surface staining:** Surface AbSeq/SMK bind to surface proteins
- 2 **Preservation:** Preserves cell status
- 3 **Permeabilization:** Allows entry of AbSeq molecule into cell


IC CITE-seq Assay chemistry

SAFE STOPPING POINT:


User will have an option to store the cells for up to 24 h in BD® OMICS-Guard Buffer before moving to permeabilization




Surface AbSeq Staining



Preservation



Permeabilization




Blocking

- 1 **Surface staining:** Surface AbSeq/SMK bind to surface proteins
- 2 **Preservation:** Preserves cell status
- 3 **Permeabilization:** Allows entry of AbSeq molecule into cell
- 4 **Blocking:** Limits noise from nonspecific binding of AbSeq Ab-oligo


IC CITE-seq Assay chemistry

SAFE STOPPING POINT:


User will have an option to store the cells for up to 24 h in BD® OMICS-Guard Buffer before moving to permeabilization




Surface AbSeq Staining




Preservation



Permeabilization



Blocking




IC AbSeq Staining

- 1 **Surface staining:** Surface AbSeq/SMK bind to surface proteins
- 2 **Preservation:** Preserves cell status
- 3 **Permeabilization:** Allows entry of AbSeq molecule into cell
- 4 **Blocking:** Limits noise from nonspecific binding of AbSeq Ab-oligo
- 5 **IC staining:** IC AbSeq Ab-oligos bind to intracellular proteins

Load cells and beads

IC CITE-seq Assay chemistry


Pair ONE cell with ONE barcoded bead in microwell



Lyse cells

IC CITE-seq Assay chemistry

Lyse cell to hybridize mRNA onto barcoded capture oligos on bead



Barcoded bead

Lysed cell

mRNA 


BD® AbSeq Ab-Oligo 

Capture oligo 

Retrieve beads




IC CITE-seq Assay chemistry




Barcoded bead
AbSeq oligos and mRNAs

Synthesize cDNA

IC CITE-seq Assay chemistry



Library preparation, sequencing and data analysis




IC CITE-seq Assay protocols

IC CITE-seq Assay



Use one of below

- Single-Cell Labeling with BD[®] AbSeq Ab-Oligos (1 plex to 40 plex)
- Single-Cell Labeling with BD[®] AbSeq Ab-Oligos (41 plex to 100 plex)
- Single-Cell Labeling with BD[®] Single-Cell Multiplexing Kit and BD[®] AbSeq Ab-Oligos (1 plex to 40 plex)
- Single-Cell Labeling with BD[®] Single-Cell Multiplexing Kit and BD[®] AbSeq Ab-Oligos (41 plex to 100 plex)



Modifications to single-cell capture workflow

Use one of below

- Single-Cell Capture and cDNA Synthesis with BD Rhapsody™ Single-Cell Analysis System 23-22951(02)
- Single-Cell Capture and cDNA Synthesis with BD Rhapsody™ Express Single-Cell Analysis System 23-22952(02)
- BD Rhapsody™ HT Single-Cell Analysis System Single-Cell Capture and cDNA Synthesis 23-24252(01)
- BD Rhapsody™ HT Xpress System Single-Cell Capture and cDNA Synthesis 23-24253(01)

Use one of below

- mRNA Whole Transcriptome Analysis (WTA) and BD[®] AbSeq Library Preparation Protocol 23-24118(02)
- mRNA Whole Transcriptome Analysis (WTA), BD[®] AbSeq, and Sample Tag Library Preparation Protocol 23-24120(02)

Sequencing and bioinformatics analysis

- BD Rhapsody™ Sequence Analysis Pipeline
- BD Rhapsody™ Single-Cell Multiomics Analysis Setup User Guide 23-21333(06)
- BD Rhapsody™ Single-Cell Multiomics Bioinformatics Handbook 23-21713(05)

IC CITE-seq Assay materials for each step

IC CITE-seq Assay

Surface
staining

Materials needed:

- Surface BD[®] AbSeq Ab-Oligos(MTS AbSeq or custom AbSeq)
- (Optional) BD[®] AbSeq IDP
- (Optional) BD[®] Human SMK if sample multiplexing



40 min



120 min

IC CITE-seq Assay materials for each step

IC CITE-seq Assay



Materials needed:

- BD Rhapsody™ Intracellular AbSeq Buffer Kit
- BD® RNase Inhibitor OR RNase Inhibitor (New England Bio, M0314L)
- BD® OMICS-Guard Sample Preservation Buffer Kit
- BD® AbSeq Enhancer Kit
- Intracellular BD® AbSeq Ab-Oligos
- (Optional) BD® Human SMK if sample multiplexing




40 min



130 min


Critical step:

While slowly vortexing the cells, add drops of ice-cold permeabilization reagent equaling 1 mL



IC CITE-seq Assay materials for each step

IC CITE-seq Assay



Materials needed:

- BD Rhapsody™ HT Xpress System Package or BD Rhapsody™ Express Single-Cell Analysis System Package
- BD Rhapsody™ Enhanced Cartridge Reagent Kit
- BD Rhapsody™ Cartridge Kit or 8-Lane Cartridge
- BD Rhapsody™ cDNA Kit
- **Proteinase K** from BD Rhapsody™ Intracellular AbSeq Buffer Kit
- (Optional) BD Rhapsody™ Scanner
- (Optional) **Vybrant™ DyeCycle™ Green** if using scanner



80 min



90 min

IC CITE-seq Assay materials for each step

IC CITE-seq Assay



Materials needed:

- BD Rhapsody™ WTA Amplification Kit

Sequencing recommendation

	ILLUMINA™ Sequencers (RUO)	Element Biosciences AVITI™ Sequencer (RUO)
Read length	Read 1: 51 cycles i7: 8 cycles i5: 8 cycles Read 2: 71 cycles	Read 1: 52 cycles i7: 8 cycles i5: 8 cycles Read 2: 72 cycles
PhiX	1%	1%
Sequencing depth	600–1,000 reads per protein specificity per cell	
Notes	<ul style="list-style-type: none"> To determine the ratio of BD Rhapsody™ WTA Assay library to BD® AbSeq Assay library to pool for sequencing, contact your local field application specialist (FAS) or scomix@bd.com to access the sequencing calculator. Avoid pooling >70% BD® AbSeq Assay library as it may impact mRNA sequence quality. 	

IC BD[®] AbSeq Ab-Oligos data interpretation

Analyze BD[®] AbSeq Ab-Oligo data using FlowJo[™] Software

BD Rhapsody[™] Sequence Analysis Pipeline output files


Output files

Output	File	Content
Metrics summary	<sample_name>_Metrics_Summary.csv	Report containing sequencing, molecules, and cell metrics
Pipeline report HTML	<sample_name>_Pipeline_Report.html	Summary report containing the results from the sequencing analysis pipeline run
BAM and BAM Index	<sample_name>.BAM <sample_name>.BAM.bai	Alignment file of R2 and associated R1 annotations
Data tables ^a	<sample_name>_RSEC_MolsPerCell_MEX.zip <sample_name>_DBEC_MolsPerCell_MEX.zip	Molecules per bioproduct per cell, based on RSEC or DBEC
	<sample_name>_RSEC_MolsPerCell_Unfiltered_MEX.zip	Unfiltered tables containing all cell labels with ≥10 reads
Bioproduct Statistics	<sample_name>_Bioproduct_Stats.csv	Metrics from RSEC and DBEC Unique Molecular Identifier adjustment algorithms on a per-bioproduct basis
Single-cell analysis tool Inputs	<sample_name>_Seurat.rds	Seurat (.rds) input file containing RSEC molecules data table and all cell annotation metadata.
	<sample_name>.h5ad	Scanpy input file containing RSEC molecules data table and all cell annotation metadata.

^a. For a multiplexed samples run, the tables contain counts for putative cells from all samples combined.

Outputs when multiplex option selected

Output	File	Content
Sample Tags metrics	<sample_name>_Sample_Tag_Metrics.csv	Metrics from the sample determination algorithm
Sample Tag calls	<sample_name>_Sample_Tag_Calls.csv	Assigned Sample Tag for each putative cell
Per-sample folder	<sample_name>_Sample_Tag<number>.zip <sample_name>_Multiplet_and_Undetermined.zip	Data tables, metric summary, and expression matrix for a particular sample. Note: For putative cells that could not be assigned a specific Sample Tag, a Multiplet_and_Undetermined.zip file is also output.





For multiplexed samples, sample tag specific MolsPerCell files can be found in the per-sample ZIP file


Contact local bioinformatics specialist or scomix@bd.com if you have any questions about file conversion

Make sure the treated and untreated samples have the same headers



FlowJo™ Software analysis



- Open a new FlowJo Workspace and drag in .csv files
- Double-click on sample to open graph window
- Click Transform Data Display drop down
 - Select “Custom Axis ...”





FlowJo™ Software analysis—Changing visualization settings (default initial starting point for all parameters)




- In Transform Window, select all parameters in “Apply to Parameters ...” list
- Change “Scale” to Biex
- Default Initial Settings (adjust as needed for individual data sets)
 - Extra Neg. Decades 0
 - Width Basis -10
 - Positive Decades 3.0
 - Max 2000
- Click Apply

FlowJo™ Software analysis—changing visualization settings (adjusted for individual parameters)


- Look at all parameters of interest in the graph window
- If data is off scale, adjust dynamic range of parameter display
 - Transform Data Display -> Customize Axis ...
 - Change Max so data is within the displayed range
 - Click Apply



Max changed from 2,000 to 20,000




Max also adjusted for y-axis




FlowJo™ Software analysis—Phenotype gating

- Gate on cells of interest
 - Example: Lymphocytes, CD3+ subset
 - PBMC sample
 - Gate out CD11c+ Monocytes
 - Gate on CD3+ Lymphocytes

Histogram plots – Interval gates



Bivariate contour plots – Rectangle gates




FlowJo™ Software analysis—IC AbSeq visualization

Bivariate plots


Histogram plot

T-bet vs CD3
Quadrant gate




PBMC-Lymph population

Granzyme B vs CD8
Quadrant gate




PBMC-L > CD3+ population

Active Caspase-3
Interval gate




PBMC-L > CD3+ population



Sequencing data analysis using the BD Rhapsody™ Sequence Analysis Pipeline

- Get access to the BD Rhapsody™ Sequence Analysis Pipeline on the [Seven Bridges Genomics Platform](#) or on a local installation
- Acquire the AbSeq reference file (.fasta) from the [BD AbSeq Panel Generator](#)
- Set up analysis following the *BD® Single-Cell Multiomics Analysis Setup User Guide (23-21333)*

Supporting you with your single-cell experiments



Getting help from single-cell experts

Visit us at scomix.bd.com to view our resource library, learning center and FAQs



In need of technical support

BD technical service support is here to help with instrument support. Contact us **email** at scomix@bd.com or online at <https://scomix.bd.com/hc/en-us/requests/new> to submit a ticket



Ordering BD[®] AbSeq Ab-Oligos and intracellular CITE-seq products

To request a quote or place an order, visit bdbiosciences.com/scM-reagents, email scomix@bd.com or contact your local BD sales representative.

Thank you



For Research Use Only. Not for use in diagnostic or therapeutic procedures.

BD, the BD Logo, BD Rhapsody and FlowJo are trademarks of Becton, Dickinson and Company or its affiliates. All other trademarks are the property of their respective owners. © 2023 BD. All rights reserved. BD-111679 (v1.0) 1223