

BD Rhapsody™ VDJ CDR3 Protocols

Overview of performance data




BD Rhapsody™ VDJ CDR3 Protocol(s)


- 1. Rhapsody VDJ CDR3 protocol** – Utilizes BD Targeted kit + Immune response panel and enables customers to run VDJ assay
- 2. Rhapsody VDJ CDR3 + BD® AbSeq protocol** - Utilizes BD Targeted kit + Immune response panel + AbSeq and enables customers to run VDJ + AbSeq
- 3. Rhapsody VDJ CDR3 + SMK protocol** - Utilizes BD Targeted kit + Immune response panel + SMK and enables customers to run VDJ on multiple samples together
- 4. Rhapsody VDJ CDR3 + BD AbSeq + SMK protocol** - Utilizes BD Targeted kit + Immune response panel + AbSeq + SMK and enables customers to run VDJ + AbSeq on multiple samples together



Overall Workflow



Workflow




New synthetic poly-A tail is added on cDNA 3' end via template switching with a poly-T TSO (template switching oligo)

Synthetic poly-A tail on cDNA 3' end hybridizes to poly-T on its own bead.

Addition of DNA polymerase and reagents allows copying of new barcodes to cDNA 3' end (equivalent to mRNA 5' end)

Ability to amplify from the 5' end of cDNA

Workflow (cont.)



Technical Features & Performance Summary

Features	Metric	Current performance
Library quality	% reads assigned to cell labels	>80 %**
	% reads aligned to amplicons	
Sequencing depth	# of reads per cell for V(D)J assembly	<5000 3000 (optimal) 2000 (met specs)
Cells processed	# of cells that can be processed	1000 - 20000
Assay time	Time of Assay (only V(D)J library)	2 working days
Workflow	1 st Stopping Point	1-3 hours
Cost, usability, multiomics capability	AbSeq / SMK	Compatible
Performance	% recovery of cells loaded	>85%**
Species	Species	Human & Mouse
Compatibility	AbSeq (40-plex), SMK Not compatible with WTA	Compatible (optimized)
Samples tested	PBMCs, enriched samples	Optimized

** - Average performance data based on limited number of internal experiments



BD Rhapsody™ VDJ CDR3 Protocol

Provided and needed components



- 1) Tested, optimized VDJ CDR3 protocols
- 2) Primer **sequences** for human CDR3 region of TCR and BCR
- 3) Bioinformatics pipeline



- 1) The Rhapsody Targeted mRNA and AbSeq reagent kit (**#633731, #633732, #633733, #633734**) along with immune response panel (**#633750**)
- 2) AbSeq or SMK (as required)
- 3) Need to order the **primers required for TCR and BCR** (listed explicitly in the BD protocols),
- 4) Order **additional reagents** required (listed explicitly in the BD protocols).



Compatibility with Targeted Assays

- While the VDJ protocol can technically work with any Targeted panel, the immune response panel with genes specific for immune cells greatly aids in immune cell calling, a feature critical for VDJ output.
- The T cell panel can work if the BD VDJ assay is being used for T cell analyses only.
- The Onco-BC panel will not be effective in cell calling for VDJ
- To work effectively with our CDR3 VDJ assay, custom panels are advised to include immune response panel markers.

BD Rhapsody™ VDJ CDR3 Protocol

Required Materials and Primer Sequences

Required and recommended materials

Required reagents

Store the reagents at the storage temperature specified on the label.

Material	Supplier	Catalog No.
BD Rhapsody™ Targeted mRNA and AbSeq Amplification Kit	BD Biosciences	633774
Agencourt® AMPure® XP magnetic beads	Beckman Coulter	A63880
100% ethyl alcohol	Major supplier	–
10 mM Tris-HCl with 0.05% Tween-20, pH 8.0 (Tris-Tween20)	Teknova	T1485
Nuclease-free water	Major Supplier	–
Template switch oligo, TSO (5' T TTT TTT TTT TTT TTT TTT TTT TTT rG rG 3'). Stock concentration of TSO is 100 μM in 10 mM Tris, 0.1 mM EDTA, pH 8.0	Major Supplier	
Klenow Fragment (3' → 5' exo-) (includes NEBuffer 2)	New England Biolabs	M0212L
10 mM dNTP	New England Biolabs	N0447L
1 M MgCl ₂	Major Supplier	–
1 M Tris-HCl, pH 8.0 (diluted to 50 mM Tris-HCl for the assay)	Major Supplier	–

Human T cell PCR2 primers

Primer Name	Primer Sequence (5' – 3')
TRAC_N2	CAGACGGTGTGCTCTCCGATCTATCAAAATCGGTGAATAGGCAGAC
TRBC_N2	CAGACGGTGTGCTCTCCGATCTGATCTCTGCTCTCTGATGGCTCA
TRDC_N2	CAGACGGTGTGCTCTCCGATCTATATCTCTGGGGTAGAATCTCTTC
TRGC_N2	CAGACGGTGTGCTCTCCGATCTGGGAAACATCTGCATCAAGTTG

Pooling Human T cell PCR2 primer panel

Name	Primer Stock Conc. (μM)	Volume per primer (μL)	DNA Suspension Buffer (μL)	Final volume (μL)
TRAC_N2 TRBC_N2 TRDC_N2 TRGC_N2	100	25.00	400.00	500

Human B cell PCR1 primers

Primer Name	Primer Sequence (5' – 3')
IGHA_N1	CACAGTCAGATCTCTGGCT
IGHD_N1	GATCTCTCTCTACTCTTGTCTGG
IGHB_N1	CGCTGAAGGTTTGTGTGTCG
IGHG_N1	TGTTGCTGGGCTTGTGAT
IGHM_N1	CGTCTCTCTCTTGTGTGCGT
IGKC_N1	TTTGTGTCTCTCTGTAGTCTGCT
IGLC_N1	TGTAGCTCTCTGTGGGACTTC

Pooling Human B cell PCR1 primer panel

Name	Primer Stock Conc. (μM)	Volume per primer (μL)	DNA Suspension Buffer (μL)	Final volume (μL)
IGHA_N1 IGHD_N1 IGHB_N1 IGHG_N1 IGHM_N1 IGKC_N1 IGLC_N1	100	15.00	420.00	525




Sensitivity and Specificity Data

Examples of TCR and BCR sensitivity and specificity metrics

**Combined data from multiple experiments and 5 donors. Values displayed are mean +/- standard deviation.




BCR Performance – Sensitivity




BCR Performance – Specificity


(% non-B cells with BCR CDR3 chain sequences)




BCR Performance – Specificity (% non-B cells with BCR CDR3 chain sequences)




TCR Performance – Sensitivity



TCR Performance – Sensitivity




TCR Performance – Sensitivity



TCR Performance – Specificity


(% non-T cells with TCR CDR3 chain sequences)




Sequencing Depth Data

An overview of different immune cell subset detection at varying sequencing depths


B cells




CD4+ T cells



CD8+ T cells



Gamma Delta T cells



BD VDJ Bioinformatics Pipeline Output

VDJ Pipeline

All reads
put into
the same
input



DRAFT BD Rhapsody™ Analysis Pipeline run - 10-16-19 01:55:29

Last update by liukai1029 on Oct. 15, 2019 18:55

App: BD Rhapsody™ Analysis Pipeline - Revision: 0

Task Inputs Execution Settings

Inputs

Batching Off

▼ **AbSeq Reference** Select file(s)

No files selected

▼ **Reads *** Select file(s)

No files selected

This field is required and cannot be empty.

▼ **Reference *** Select file(s)

No files selected

This field is required and cannot be empty.

App Settings

Edit parameters Show editable ▾

▶ **Putative_Cell_Calling_Settings** (#Putative_Cell_Calling_Settings)

▶ **Subsample_Settings** (#Subsample_Settings)

▼ **VDJ_Settings** (#VDJ_Settings)

VDJ Species Version

- ✓ No value
- Human VDJ - BCR and TCR
- Human VDJ - BCR only
- Human VDJ - TCR only
- Mouse VDJ - BCR and TCR**
- Mouse VDJ - BCR only
- Mouse VDJ - TCR only

Subsample Sample Tags

No value



_Metrics Summary File

#Sequencing Quality#						
Total_Reads_in_FASTQ	Pct_Reads_Too_Short	Pct_Reads_Low_Base_Quality	Pct_Reads_High_SNF	Pct_Reads_Filtered_Out	Total_Reads_After_Quality_Filtering	Library
48068574	1.1	6.98	6.08	11.2	42683736	BCR10k
61090421	0.71	5.81	3.55	9.16	55493393	TCR10k
14765915	0.11	11.51	5.66	15.59	12464008	mRNA10k
123924910	0.79	6.94	4.78	10.72	110641137	Combined_stats

#Library Quality#						
Total_Filtered_Reads	Pct_Contaminating_PhiX_Reads_in_Filtered_R2	Pct_Q30_Bases_in_Filtered_R2	Pct_Assigned_to_Cell_Labels	Pct_Cellular_Reads_Aligned_Uniquely_to_Amplicons	Library	
42683736	0.03	62.61	84.39	79.76	BCR10k	
55493393	0.07	64.58	82.17	76.58	TCR10k	
12464008	0.01	66.32	90.6	86.88	mRNA10k	
110641137	0.05	64.01	83.98	78.97	Combined_stats	

_VDJ_Metrics Summary File

#Overall VDJ Metrics#									
Reads_Cellular_Aligned_to_VDJ	Reads_CDR3_Valid_Unfiltered	Reads_CDR3_Valid_Putative	Pct_Reads_CDR3_Valid_from_Putative_Cells	Reads_CDR3_Valid_Putative_Corrected	Pct_Reads_CDR3_Valid_Corrected_from_Putative_Cells	Mean_Reads_CDR3_Valid_Corrected_per_Putative_Cell	Molecules_Unfiltered	Molecules_Corrected_Putative	Mean_Molecules_Corrected_per_Putative_Cell
63975808	48826575	38366518	78.58	36652359	75.07	45082.85	508537	23366	28.74

#Chain Type Metrics#									
Chain_Type	Reads_CDR3_Valid_Unfiltered	Reads_CDR3_Valid_Putative	Pct_Reads_CDR3_Valid_from_Putative_Cells	Reads_CDR3_Valid_Putative_Corrected	Pct_Reads_CDR3_Valid_Corrected_from_Putative_Cells	Mean_Reads_CDR3_Valid_Corrected_per_Putative_Cell	Molecules_Unfiltered	Molecules_Corrected_Putative	Mean_Molecules_Corrected_per_Putative_Cell
BCR_Heavy	15201582	12195119	80.22	11850057	77.95	14575.72	120649	6220	7.65
BCR_Kappa	4046534	2970188	73.4	2701916	66.77	3323.39	43742	2576	3.17
BCR_Lambda	4823527	4013658	83.21	3948324	81.86	4856.49	44895	2999	3.69
TCR_Alpha	11672066	8742143	74.9	8004555	68.58	9845.7	132640	4406	5.42
TCR_Beta	9801908	7541670	76.94	7420795	75.71	9127.67	125525	5588	6.87
TCR_Delta	808033	657997	81.43	657120	81.32	808.27	11159	356	0.44
TCR_Gamma	2472925	2245743	90.81	2069592	83.69	2545.62	29927	1221	1.5

#Cell Type Metrics#									
Cell_Type_Experimental	Number_cells	BCR_Paired_Chains_Percent	TCR_Paired_Chains_Percent	BCR_Heavy_Percent_Cells_Positive	BCR_Heavy_Mean_Molecules_per_Cell	BCR_Kappa_Percent_Cells_Positive	BCR_Kappa_Mean_Molecules_per_Cell	BCR_Lambda_Percent_Cells_Positive	BCR_Lambda_Mean_Molecules_per_Cell
B	92	94.57	0	94.57	56.2	53.26	26.98	44.57	17.09
Dendritic	12	0	0	0	0	0	0	0	0
Monocyte_classical	128	0.78	0.78	0.78	0.14	0.78	0.59	0	0
Monocyte_nonclass	17	0	0	0	0	0	0	0	0
Natural_killer	210	1.43	3.33	2.86	4.74	1.43	0.09	2.38	6.69
T_CD4_memory	183	0	84.15	0.55	0.01	0	0	0	0
T_CD4_naive	46	0	84.78	2.17	0.04	2.17	0.02	0	0
T_CD8_memory	55	1.82	81.82	3.64	0.47	0	0	1.82	0.42
T_CD8_naive	55	0	89.09	1.82	0.13	0	0	0	0
T_gamma_delta	15	0	86.67	0	0	0	0	0	0



_VDJ_perCell File

Cell_Index	Total_VDJ_Read_Count	Total_VDJ_Molecule_Count	BCR_Heavy_V_gene_Dominant	BCR_Heavy_D_gene_Dominant	BCR_Heavy_J_gene_Dominant	BCR_Heavy_C_gene_Dominant	BCR_Heavy_CDR3_Nucleotide_Dominant	BCR_Heavy_CDR3_Translation_Dominant	BCR_Heavy_Read_Count	BCR_Heavy_Molecule_Count
195260	46123	96	IGHV4-34*01	IGHD6-19*01	IGHJ5*02	human_IGHG	GCGAGAGGAGGGA	(ARGGGYRSGWYGA	16914	27
275040	29983	82	IGHV4-39*01	IGHD3-16*01	IGHJ5*02	human_IGHM	GCGAGACGTGATCARRDRLGS		24035	58
790251	22513	87	IGHV3-30*02	IGHD3-3*01	IGHJ4*02	human_IGHM	GCGAAAGATTGGA	AKDWSAFDY	346	16
856783	4172	24							0	0
565345	15494	40	IGHV3-15*01	IGHD4-17*01	IGHJ6*02	human_IGHD	ACCACAGAGCTCGTTELDYGVGGYYYY		11283	24
512464	0	0							0	0
251681	0	0							0	0
558698	22957	59	IGHV4-4*06	IGHD1-7*01	IGHJ5*02	human_IGHM	GCGAGAGTTCATA	ARVHNWNYPADDN	18124	45
11337	0	0							0	0
777579	26444	83	IGHV4-31*02	IGHD2-2*01	IGHJ3*02	human_IGHM	GCGAGGTATTGTA	ARYCSSPSCSGAFD	18452	38
734829	44185	121	IGHV1-2*02	IGHD6-19*01	IGHJ4*02	human_IGHM	GCGAGTGGTGGC	ASGGWYHY	33403	91
853307	18633	70	IGHV5-51*01	IGHD4-17*01	IGHJ3*02	human_IGHA	GCGAGACTCTCGT	ARLSSPDYDLPLGR	13215	43
831417	0	0							0	0
6245	0	0							0	0
695148	31532	75	IGHV3-30*01	IGHD6-19*01	IGHJ4*02	human_IGHM	GCGAGATTGTATA	ARLYSSGWGHVDY	13888	43
688461	0	0							0	0
673277	51375	94	IGHV1-8*01	IGHD3-16*01	IGHJ5*02	human_IGHM	GCGAGAGGCGCC	(ARGALLSGAHWGY	40231	63
813557	27089	62	IGHV1-18*01	IGHD5-12*01	IGHJ4*02	human_IGHM	GCGAGACCGGGG	(ARPGGDYSGYDYG	21345	39
592931	25211	71	IGHV1-2*02	IGHD2-2*02	IGHJ6*02	human_IGHG	GCGAGGGAAGCTC	AREAVVLPAAIRHFC	7106	15



_VDJ_perCell File

Cell_Index	Total_VDJ_Read_Count	Total_VDJ_Molecule_Count	BCR_Heavy_V_gene_Dominant	BCR_Heavy_D_gene_Dominant	BCR_Heavy_J_gene_Dominant	BCR_Heavy_C_gene_Dominant	Cell_Index	BCR_Paired_Chains	TCR_Paired_Chains	Cell_Type_Experimental	Heavy_Molecule_Count
195260	46123	96	IGHV4-34*01	IGHD6-19*01	IGHJ5*02	human_IGHG	195260	TRUE	FALSE	B	27
275040	29983	82	IGHV4-39*01	IGHD3-16*01	IGHJ5*02	human_IGHM	275040	TRUE	FALSE	B	58
790251	22513	87	IGHV3-30*02	IGHD3-3*01	IGHJ4*02	human_IGHM	790251	TRUE	FALSE	B	16
856783	4172	24					856783	FALSE	TRUE	T_CD8_memory	0
565345	15494	40	IGHV3-15*01	IGHD4-17*01	IGHJ6*02	human_IGHD	565345	TRUE	FALSE	B	24
512464	0	0					512464	FALSE	FALSE	Natural_killer	0
251681	0	0					251681	FALSE	FALSE	Natural_killer	0
558698	22957	59	IGHV4-4*06	IGHD1-7*01	IGHJ5*02	human_IGHM	558698	TRUE	FALSE	B	45
11337	0	0					11337	FALSE	FALSE	Natural_killer	0
777579	26444	83	IGHV4-31*02	IGHD2-2*01	IGHJ3*02	human_IGHM	777579	TRUE	FALSE	B	38
734829	44185	121	IGHV1-2*02	IGHD6-19*01	IGHJ4*02	human_IGHM	734829	TRUE	FALSE	B	91
853307	18633	70	IGHV5-51*01	IGHD4-17*01	IGHJ3*02	human_IGHA	853307	TRUE	FALSE	B	43
831417	0	0					831417	FALSE	FALSE	Monocyte_classical	0
6245	0	0					6245	FALSE	FALSE	Monocyte_classical	0
695148	31532	75	IGHV3-30*01	IGHD6-19*01	IGHJ4*02	human_IGHM	695148	TRUE	FALSE	B	43
688461	0	0					688461	FALSE	FALSE	Natural_killer	0
673277	51375	94	IGHV1-8*01	IGHD3-16*01	IGHJ5*02	human_IGHM	673277	TRUE	FALSE	B	63
813557	27089	62	IGHV1-18*01	IGHD5-12*01	IGHJ4*02	human_IGHM	813557	TRUE	FALSE	B	39
592931	25211	71	IGHV1-2*02	IGHD2-2*02	IGHJ6*02	human_IGHG	592931	TRUE	FALSE	B	15

VDJ + SMK

An overview of VDJ + SMK performance data performed using 2 donor samples

Sample Tag Metrics




Sensitivity – 99.33%

Specificity



B	99.16
CD4T_memory	99.76
CD4T_naive	99.76
CD8T_memory	100
CD8T_naive	98.85
DC	100
Monocyte_classical	100
Monocyte_nonclassical	100
NK	97.95
gdT	100

BCR Performance – Sensitivity




BCR Performance – Specificity

(% non-B cells with BCR CDR3 chain sequences)




BCR Performance – Specificity


(% non-B cells with BCR CDR3 chain sequences)




TCR Performance – Sensitivity



TCR Performance – Sensitivity




TCR Performance – Sensitivity



TCR Performance – Specificity









(% non-T cells with TCR CDR3 chain sequences)




VDJ + BD[®] AbSeq

An overview of VDJ + AbSeq performance data


BD AbSeq Marker Detection in B Cells



BD AbSeq Marker Detection in T Cells





BCR Performance – Sensitivity




BCR Performance – Specificity

(% non-B cells with BCR CDR3 chain sequences)





BCR Performance – Specificity


(% non-B cells with BCR CDR3 chain sequences)



TCR Performance – Sensitivity




TCR Performance – Sensitivity



TCR Performance – Specificity

(% non-T cells with TCR CDR3 chain sequences)



Mouse Primer Panel

TCR Primer Panel

Mouse T cell PCR1 Primers

Primer Names	IDT Primer Stock Conc. (μM)	Vol/Primer Required (μl)	Total Vol of IDTE (μl)	Final Vol (μl)	Sequence
TRAC N1	200	12.5	437.5	500	TTTTCGGCACATTGATTTGGGAG
TRBC N1					CTCAGGCAGTAGCTATAATTGCT
TRDC N1					CAATCTTCTTGGATGATCTGAGACT
TRGC1-TRGC2 N1					GGAAAGAACTTTTCAAGGAGACAAAGG
TRGC4 N1					CATCCTTTTCTTTCCAATACACCC

Mouse T cell PCR2 Primers

TRAC N2	200	12.5	437.5	500	CAGACGTGTGCTCTTCCGATCTAGGTTCTGGGTTCTGGATGT
TRBC N2					CAGACGTGTGCTCTTCCGATCTCAATCTCTGCTTTTGATGGCTC
TRDC N2					CAGACGTGTGCTCTTCCGATCTGTAGAAATCTTTCACCAGACAAGC
TRGC1-TRGC2 N2					CAGACGTGTGCTCTTCCGATCTTTGGGGGAAATGTCTGCA
TRGC4 N2					CAGACGTGTGCTCTTCCGATCTAATAGTAGGCTTGGGAGAAAAGTCTG

BCR Primer Panel

Mouse B cell PCR1 Primers

Primer Names	IDT Primer Stock Conc. (μM)	Vol/Primer Required (μl)	Total Vol of IDTE (μl)	Final Vol (μl)	Sequence
IGHA N1	200	7.5	442.5	525	AACTGGCTGCTCATGGTGTA
IGHD N1					AAGTGTGGTTGAGGTTTCAGTTCTG
IGHE N1					GAAGTTCACAGTGCTCATGTTC
IGHG1 N1					CAGAGTGTAGAGGTCAGACT
IGHG2A-IGHG2C N1					TCGAGGTTACAGTCACTGAG
IGHG2B N1					GATCCAGAGTTCCAAGTCACAG
IGHG3 N1					TACGTTGCAGATGACAGTCT
IGHM N1					TGGATGACTTCAGTGTGTTCTG
IGKC N1					TGTAGGTGCTGTCTTTGCTG
IGLC1 N1					CTGTAAGTCTATGCCTTTCCC
IGLC2-IGLC3 N1					TTGGTGGGATTTGAAGTGTCC

Mouse B cell PCR2 Primers

IGHA N2	200	39.0	91.0	520	CAGACGTGTGCTCTCCGATCTTGTCAGTGGGTAGATGGTGG
IGHD N2					CAGACGTGTGCTCTCCGATCTCTGACTTCCAATTAACAGCC
IGHE N2					CAGACGTGTGCTCTCCGATCTTAGAGCTGAGGGTTCCTGATAG
IGHG1 N2					CAGACGTGTGCTCTCCGATCTCAGTGGATAGACAGATGGGGGT
IGHG2A-IGHG2C N2					CAGACGTGTGCTCTCCGATCTATGGGGCTGTTGTTTTGG
IGHG2B N2					CAGACGTGTGCTCTCCGATCTGTGGATAGACTGATGGGGGTGTT
IGHG3 N2					CAGACGTGTGCTCTCCGATCTAGGGAAGTAGCCTTTGACAAG
IGHM N2					CAGACGTGTGCTCTCCGATCTGACATTTGGGAAGGACTGACTC
IGKC N2					CAGACGTGTGCTCTCCGATCTAGATGTTAACTGCTCACTGGATG
IGLC1 N2					CAGACGTGTGCTCTCCGATCTGTTAGTCTCGAGCTCTTCAGA
IGLC2-IGLC3 N2					CAGACGTGTGCTCTCCGATCTCAGTGTGGCTTTGTTTTCT





For research only. Not for use in diagnostic or therapeutic procedures.

23-22449-00

BD, the BD Logo, Rhapsody and SeqGeq are trademarks of Becton, Dickinson and Company or its affiliates. All other trademarks are the property of their respective owners. © 2020 BD. All rights reserved.