

BD Rhapsody™ Cartridge Reagent Kit

BD Biosciences Kit Part Number: 633731

Purpose

The purpose of this kit is to provide the necessary reagents to support library preparation.

Contents

PN	Component	Cap color	Quantity	Volume per unit	Storage
650000060	Cartridge Wash Buffer 1	Neutral	1 bottle	7 mL	2°C to 8°C
650000061	Cartridge Wash Buffer 2	Neutral	1 bottle	4 mL	
650000062	Sample Buffer	Neutral	1 bottle	28 mL	
650000064	Lysis Buffer	Neutral	4 bottles	15 mL	
650000065	Bead Wash Buffer	Neutral	1 bottle	10 mL	
650000063	1 M DTT	White	1 vial	400 µL	
650000089	Cell Capture Beads	Brown	4 vials	2 mL	
650000090	Waste Collection Container	Neutral	4 each	—	

Handling

Store the reagents at the specified storage temperature. Before using reagents, see the *BD Rhapsody™ Single-Cell Analysis System User Guide* (Doc ID: 47395).

Safety

Use standard laboratory safety protocols. Read and understand the safety data sheets (SDSs) before handling chemicals. To obtain SDSs for chemicals ordered from BD Biosciences, go to regdocs.bd.com or contact BD Biosciences technical support at researchapplications@bd.com.

BD Biosciences

1.877.232.8995, prompt 2, 2 | researchapplications@bd.com | 2350 Qume Drive, San Jose, CA 95131 USA

Conditions: The information disclosed herein is not be construed as a recommendation to use the above product in violation of any patents. BD Biosciences will not be held responsible for patent infringement or other violations that may occur with the use of our products. Purchase does not include or carry any right to resell or transfer this product either as a stand-alone product or as a component of another product. Any use of this product other than the permitted use without the express written authorization of Becton, Dickinson and Company is strictly prohibited.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Not for resale.

© 2018 BD. BD and the BD Logo are trademarks of Becton, Dickinson and Company. 04/2018

BD Rhapsody Cartridge Reagent Kit Product Insert

Doc ID: 54141 Rev. 2.0

